Product data sheet

wiroSPADE

WOODEN TONGUE DEPRESSOR

Smoothed | Rounded ends

PREMIUM+



Product data sheet

wiroSPADE

Manufacturer

Company + Address	WIROS Wilfried Rosbach GmbH, Karl-Arnold-Str. 5, 47877 Willich, Germany
SRN	DE-MF-000004956
Trademark	care & serve®

Authorisation

CH REP	Swiss AR Services AG, Industriestrasse 47, 6300 Zug, Switzerland

Product description

Category	PREMIUM+
Product description	care & serve® PREMIUM+ | wiroSPADE | Wooden tongue depressor Smoothed | Rounded ends
Product features	Smoothed Rounded ends

Material

Material	Birchwood

Application

Intended use	Disposable wooden medical instrument to depress the tongue and expand the cheeks to get a better view of the oral cavity and throat.
Scope of application	Pushing down the tongue to inspect the oral cavity Folding back the cheeks to inspect the mucous membranes etc.
Application group	Use is restricted to appropriately trained professionals.
Patient group	Wooden tongue depressors suitable for all patient groups.
Indication	Pushing down the tongue to inspect the oral cavity Folding back the cheeks to inspect the mucous membranes etc.
Contraindication	Not suitable for use when treating immunocompromised patients or premature infants. Non-sterile wood products should not be used for the collection, handling and processing of samples for microbiological testing, wound care or in intensive care units.
Notes on use	Disposable product Not sterile
Limitation	Natural material – Please carry out visual inspection before use.
Note	All serious incidents directly related to this product must be reported to the manufacturer and the competent authorities.
Application areas	Hospital | Old people's home | inc. specialist circles Medical practices Laboratory Cosmetics | Wellness | Beauty industry
Supplementary material info	Latex-free

Product data sheet

wiroSPADE

Standards

CE marking	Class I medical device in accordance with EU 2017/745
SRN	DE-MF-000004956
Basis UDI DI	4051642-wiroSPADE-EU
	CE marking The CE mark guarantees the free movement of products and goods within the European Economic Community. All products with the CE mark comply with the essential requirements of the European regulations specified.
VO (EU) 2023/988	Regulation (EU) No 2023/988 on general product safety

Use

Instructions for use	Before use Please ensure that the disposable tongue depressor has no chips or broken areas. Only use intact and unused disposable tongue depressors. Use Use the disposable tongue depressor to carefully fix and press on the patient's tongue. During use, make sure that the depressor is not inserted too far as doing so could trigger a gag reflex. The disposable tongue depressor is then disposed of.
Safety notice	Natural material – check for chips and broken areas before use.
Disposal	Unused and uncontaminated products can be incinerated in an environmentally friendly manner or disposed of at a waste disposal site. Contaminated products must be disposed of in accordance with the national laws and regulations valid for the respective type of contamination.
Shelf life	5 years from the production date, if unused and stored correctly.
Storage	Keep out of direct sunlight. Store in a dry place. Incorrect storage can lead to a reduction in the shelf life.
Traceability	Batch-related traceability from the determination of requirements to issuing. Product data sheet Version 2 Valid from Batch 2617 +

Packaging

Recycling	Licensed in accordance with VerpackG (German Packaging Act), Licence number DE1618662887322. Please recycle the packaging.

www.wiros.de

WEB	https://www.wiros.de/shop/medizin-pflegebedarf/mundspatel/wirospade/

Variant 1

Sales unit	5000 Unit/Cardboard box
Packaging level 1	100 Unit/Polybag
Minimum order quantity	1 Cardboard box

Item numbers

Item number	Colour	Size (W x L)	Tolerance
95800500	natural	≈ 1,7 x 15 cm	± 2 mm

GTINs

Item number	Unit of Use UDI	Packaging level 1 UDI	Packaging level 2 UDI	Cardboard box UDI
95800500	4051642003825	4051642003290	----------	4051642503295

Product data sheet

wiroSPADE

Variant 2

Sales unit	5000 Unit/Cardboard box
Packaging level 1	100 Unit/Dispenser box
Minimum order quantity	1 Cardboard box

Item numbers

Item number	Colour	Size (W x L)	Tolerance
95800600	natural	≈ 1,7 x 15 cm	± 2 mm

GTINs

Item number	Unit of Use UDI	Packaging level 1 UDI	Packaging level 2 UDI	Cardboard box UDI
95800600	4051642003832	4051642003306	----------	4051642503301

The numerical data for the measurements, weights, thicknesses and tolerances are target values, and there may be some deviation.
The information represents the current state of our knowledge. It offers no guarantee of the product characteristics and does not form the basis of a contractual legal understanding. Please reconfirm any information that is older than 2 years.
Typographical and other errors excepted.
The product data sheets provided are the property of WIROS Wilfried Rosbach GmbH. WIROS Wilfried Rosbach GmbH is the owner of all copyrights and rights of use and exploitation for the product data sheets. It is not permitted to change the product data sheets or to use the pictograms. The product data sheets provided may only be used in conjunction with original WIROS products.