Combur[5] Test HC

cobas®

REF 11896954176

10

English

Intended use

The Combur[5] Test HC are five-parameter test strips for in vitro semi-quantitative determination of glucose, leukocytes, nitrite, protein and erythrocytes in urine by visual reading. These measurements are useful in the evaluation of renal, urinary and metabolic disorders.
Combur[5] Test HC are test strips for single use only.
Combur[5] Test HC are screening tests for patients. The test results can aid physicians in the diagnosis of potential pathological conditions.
The test is intended for self-testing.

Summary

Combur[5] Test HC is not intended for a near-patient testing environment.
Combur[5] Test HC provides a fast and easy means of early detection and clinical monitoring of diabetes mellitus and diseases of the kidneys and urinary tract.
Combur[5] Test HC supports your physician's efforts to detect these diseases in their early stages. It does this by indicating the presence of substances excreted in the urine, thereby providing evidence of disease even when no symptoms are present.

- **Test your urine for 2 to 3 days in succession.**
- **Self-testing cannot, however, be used in place of testing performed by a physician.**
- **Continue testing every 3 months using 2 to 3 test strips, even if your test results are normal or your physician does not detect any abnormalities. This will provide additional security.**
- **Please contact your physician if you detect + or ++ results on one or more test parameters. If you are pregnant, persistent or greater than trace amounts of glucose in urine should be investigated by your treating physician. Please test your urine again if the test results are unclear. In case all parameters are 0 after the first measuring, please test your urine for 1 or 2 more days and show the results to your physician.**
- **There is a possibility that test results might be false-positives or false-negatives. In case of doubt, please test your urine again.**
- **Please contact your physician if you experience symptoms, typical for urinary tract infections:**
 - frequent and/or painful urination (dysuria), sudden onset of the need to urinate, urine looking cloudy, blood in urine (hematuria), and /or pelvic pain.
 - Especially, if you experience fever, chills, nausea, vomiting or back pain, which could be signs for a complicated urinary infection.
- **Always consult your physician before making any decision that could affect your treatment.**

Precautions and warnings

The stopper of the test strip vial contains a non-toxic silicate-based desiccant, which must not be removed. If ingested by accident, drink large quantities of water.
For in vitro diagnostic use for lay users.
Dispose all waste materials in accordance with local rules. Prevent release into the environment, drainage system or water bodies.
Environmental hazards:
Apply all relevant local disposal regulations to determine the safe disposal.
Safety data sheet available for professional user on request.

Storage and stability

Store the package at 2-30 °C. When stored in the original vial the test strips are stable up to the expiration date specified on the box and color label.
Tightly re-cap the vial immediately after removing a test strip.
Do not use test strips if they are past the expiration date printed on the vial and color label.
Discard expired test strip.

Specimen collection and preparation

Use only clean, well-rinsed vessels to collect urine.
Do not add preservatives to the urine.
Use fresh urine that has not been centrifuged.[1] The urine specimen should not stand for more than 2 hours before testing. For specimen collection and preparation only use suitable tubes or collection containers, as false-positive readings, particularly for glucose and protein, can result from residues of detergent or strongly oxidizing disinfectants in the specimen collection vessel.[1]
Using midstream urine is recommended to avoid contamination by commensal urethral flora in both sexes.[2] Do not expose urine specimens to sunlight as this induces oxidation of bilirubin and urobilinogen and hence leads to artificially low results for these 2 parameters.[2] Vaginal secretion or menstrual blood may contaminate urine from females.[2]
Diagnosis or therapy should never be based on one test result alone but should be established in the context of all other medical findings. In doubtful cases, it is therefore advisable to repeat the test after discontinuation of the medication.

Materials provided

- 10 test strips

Materials required (but not provided)

- A clean cup for collecting urine

Directions for use

For optimum performance of the test, follow the directions given below:

- Allowing the urine to remain in the bladder for long periods of time (4-8 hours, ideally overnight) is recommended, as a means of ensuring a high level of accuracy.
- In special cases, when measuring urinary glucose, note that, in some cases, increased glucose values can only be detected after a meal. As such, performing the test 2 hours after a meal, preferably after breakfast, is recommended. If possible, this meal should include 3-4 slices of bread or 1-2 rolls.
- ⚠ In women, this test can yield false results for blood if carried out 3 days before and up until 3 days after menses; as such, women are advised against performing this test during this period.
- People with visual or color resignation impairment e.g. red-green visual impairment may misinterpret test results and therefore cannot evaluate the test unaided.

Before you test:

1. Always test your first morning urine (except for special cases of measuring urinary glucose, see above). Wash your genitals carefully and rinse thoroughly with water (do not dry). To collect urine use only clean, well-rinsed vessels.
2. Allow some urine to flow into the toilet. Collect the rest of the stream in the specimen cup.
3. Ensure that the test is performed at a temperature range of +18 °C to +32 °C and a relative humidity of 30 % to 80 %.

Performing the test:


Fig. 1 Fig. 2 Fig. 3 Fig. 4 Fig. 5 Fig. 6

1. Remove test strip. Immediately after you take a test strip out of the vial, seal the vial with the original desiccant stopper. Failure to do this could result in false results, as moisture in the vial could potentially cause color changes in the test parameters (fig. 1).
2. Compare the dry test strip against the control colors on the label (column C). Use only test strips with comparable test pad colors (fig. 2).
3. Briefly dip the test strip into the urine (approximately 1 sec.), making certain to wet all of the test parameters (fig. 3).
4. When removing the test strip, wipe its edge along the rim of the specimen cup in order to remove any excess urine (fig. 4).
5. All test parameters should be read at 1 minute (up to 120 seconds for the leukocyte test area for not clearly assignable results). Compare the colors of the test parameters on the strip with the colors on the label (fig. 5).
 Separate sets of color blocks are given for erythrocytes and hemoglobin. Compare the 5th (blood) test parameter area with both rows of colors.
 Perform the test under sufficient natural light conditions e.g. daylight.

Result sheet:

1. If the color changes in one or more test parameters, place an x in the box next to the line on the results sheet (0, + or ++) that best corresponds to the color of the corresponding test parameter (fig. 6). If a test parameter appears darker than ++, check the box under ++.
2. Enter date and your name into results sheet and take it with you to your physician.

Any color changes appearing only along the edges of the test areas, or developing after more than 2 minutes, do not have any diagnostic significance.

Additional information for Healthcare Professionals

Reagents

Each test contains per 1 cm² reactive paper area the following:
Glucose: 3,3',5,5'-tetramethylbenzidine 103.5 µg; GOD 6 U, POD 35 U
Leukocytes: Indoxylcarbonic acid ester 15.5 µg; methoxymorpholinobenzene diazonium salt 5.5 µg
Nitrite: 3-hydroxy-1,2,3,4-tetrahydro-7,8-benzoquinoline 33.5 µg; sulfanilamide 29.1 µg
Protein: 3',3'',5',5''-tetrachlorophenol-3,4,5,6-tetrabromosulfophthalein 13.9 µg
Blood: 3,3',5,5'-tetramethylbenzidine 52.8 µg; 2,5-dimethyl-2,5-dihydroperoxyhexane 297.2 µg

Test principle

Test parameters on each strip contain indicator reagents. As soon as the test parameters have been wetted with urine, a chemical reaction takes place that causes the color of the test parameters to change in relation to the concentration of the substance in question. The color change can then be assessed by comparing it with the color scale given on the vial label.
Glucose (GLU): The glucose determination is based on the specific glucose-oxidase/peroxidase reaction (GOD/POD method).
Leukocytes (LEU): The test reveals the presence of granulocyte esterases. These esterases cleave an indoxyl ester, and the indoxyl so liberated reacts with a diazonium salt to produce a violet dye.
Nitrite (NIT): The test is based on the principle of the Griess test and is specific for nitrite. The reaction reveals the presence of nitrite and hence indirectly nitrite-forming bacteria in the urine by a pink-to-red coloration of the test parameter. Even a slight pink coloration is indicative of significant bacteriuria.
Protein (PRO): The test is based on the principle of the protein error of a pH indicator. It is particularly sensitive to albumin.
Blood (ERY/Hb): The peroxidase-like action of hemoglobin and myoglobin specifically catalyzes the oxidation of the indicator by means of the organic hydroperoxide contained in the test paper to give a blue-green coloration.

Limitations - interference

Therapeutic drugs and endogenous substances were tested for a potential interference to the test parameters of the Combur[5] Test HC. All parameters were tested with negative urine samples and samples spiked to the first positive concentration range. Therapeutic drugs were tested at concentrations in urine occurring under medication with the therapeutic dosage and above. There are no significant therapeutic drug interferences up to the concentrations as presented below:

Parameter	Therapeutic drug	No interference up to	Effect above stated concentration
GLU	Ascorbic acid	2000 mg/L	False negative results
LEU	Phenazopyridine	3 mg/L	False negative results
	Salicyl uric acid	100 mg/L	False negative results
	Tetracycline	100 mg/L	False negative results
NIT	Ascorbic acid	1000 mg/L	False negative results
	Phenazopyridine	5 mg/L	Not assessable results[a)]
	Salicyl uric acid	5000 mg/L	False negative results
PRO	Gabapentin	8000 mg/L	False positive results
ERY	Ascorbic acid	800 mg/L	False negative results
	Ibuprofen	1000 mg/L	False negative results

a) not assessable results: A visual determination might not be possible for negative or low positive results due to intrinsic color of the specimen.
There are no significant endogenous substance interferences up to the concentrations as presented below:

Parameter	Endogenous substance	No interference up to	Effect above stated concentration
GLU	Urea	88780 mg/L	False negative results
	Uric acid	305 mg/L	False negative results
	Urobilinogen	200 mg/L	Not assessable results[b)]
LEU	Bilirubin	102 mg/L	Not assessable results[b)]
	a-D(+)-Glucose	60038 mg/L	False negative results
	Urobilinogen	120 mg/L	Not assessable results[b)]
NIT	Bilirubin	62 mg/L	False negative results
	Urobilinogen	400 mg/L	Not assessable results[b)]
PRO	Creatinine	10599 mg/L	False positive results
	Hemoglobin	100 mg/L	False positive results and elevated positive results
	Urea	88290 mg/L	False positive results and elevated positive results
	Urobilinogen	1000 mg/L	Not assessable results[b)]
ERY	Urobilinogen	500 mg/L	False negative and not assessable results[b)]

b) not assessable results: A visual determination might not be possible for negative or low positive results due to intrinsic color of the specimen.

Common limitations

The test should not be performed within 3 days of treatment with antibiotics and/or chemotherapy. Other factors, such as intense jogging, can increase erythrocyte (blood) and protein values without indicating illness.
Leukocytes: A negative Leukocyte test result does not rule out urinary tract infections (UTI) in every case. If you had a renal transplantation or you are pregnant and have risk factors, like Diabetes mellitus, previous pyelonephritis (type of urinary tract infection that also affects the kidney), previous premature births or late pregnancy loss, please refer to your physician for screening for bacteria in the urine.
Nitrite: Prolonged urinary retention in the bladder (4-8 hours) is essential in order to obtain an accurate result.[2] Administration of antibiotics or chemical drugs should be discontinued 3 days before the test.[3] More than 80 % of all bacteria responsible for urinary tract infections are Gram-negative rods (E.coli, Klebsiella, Enterobacter and Proteus species).[4] Most gram-negative bacteria have the ability to reduce urinary nitrate to nitrite and can therefore be detected indirectly with the test strips.[2] Normal nutrition as a rule ensures a sufficiently high content of nitrate in the urine for the detection of bacteria.[5] Some common uropathogens, e.g. Enterococcus spp. and Staphylococcus spp. (5-15 % of bacteria responsible for urinary tract infections)[4], do not reduce urinary nitrate to nitrite and will therefore not be detected whatever their urinary concentration.[2] False-negative results may occur as a result of strong diuresis with frequent voiding of urine, insufficient nitrate intake or too short retention of urine in the bladder.[2]
Attention: Nitrogen oxides present in the atmosphere may have an influence on the stability of the nitrite test parameter.[5]
Protein: False-positive readings may be found after infusion of polyvinylpyrrolidone (blood substitute).
Blood/ERY: In women the test for blood may be falsified from 3 days before up until 3 days after a period. It is therefore advisable not to perform the test during this time. After physical activity, e.g. strenuous jogging, raised values for erythrocytes and protein may occur without being signs of disease.[7]
For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Expected and result values

Based on literature. Current medical guidelines are leading.

Parameter	Expected values	Additional information
GLU	< 25 mg/dL < 1.4 mmol/L[8]	For daytime urine. Using semi-quantitative reagent strips, expected values in a healthy population are negative.[9]
LEU	< 10 Leu/µL[2]	10-100 Leu/µL borderline[2]
NIT	< 1 µmol (< 0.005 mg/dL)[10]	A positive result is indicative of urinary tract infection (UTI), but a negative result does not rule out UTI.[6]
PRO	≤ 30 mg/dL[11]	> 30 mg/dL proteinuria[11]
ERY	< 18 Ery/µL (< 3 Ery/HPF)[12]	Hematuria ≥ 18 Ery/µL (≥ 3 Ery/HPF)[13,14]
	Conversion factor 5.8 to translate chamber counting high-power field (HPF) into µL[2]	

Parameter	Result values
GLU	norm., 50, 100 mg/dL norm., 2.8, 5.6 mmol/L 0, +, ++
LEU	neg., ~ 75, ~ 500 Leu/µL 0, +, ++
NIT	neg., pos., high pos. 0, +, ++
PRO	neg., 30, 100 mg/dL neg., 0.3, 1 g/L 0, +, ++
ERY	neg., ~ 5-10, ~ 50 Ery/µL 0, +, ++

In general, urine testing can be performed for all patients, where it is possible to collect a sample.[15] Stated expected values are based on literature. All other verifications of the Combur[5] Test HC were performed using adult patient samples. For testing other patient populations (e.g. children) please contact your physician.

Specific performance data

Representative performance data are given below. Results obtained in individual laboratories may differ. The values for neg. and pos. indicate the proportion of concordant negative or positive results.
The values specified for the **limit of detection** are defined as the concentration of the analyte, which leads to a positive result in ≥ 90 % of the examined urine samples.
The **method comparison** data for visual reading are based on the comparison with the instrument **cobas** u 411 with Combur[10] Test M using at least 232 clinical samples per parameter. All test parameters were covered.

Parameter	Limit of detection	Method comparison[c)]
GLU	25-45 mg/dL	neg.: 96 %, pos.: 100 %
LEU	25-75 Leu/µL	neg.: 100 %, pos.: 62 %
NIT	0.03-0.09 mg/dL	neg.: 95 %, pos.: 91 %
PRO	10-18 mg/dL (albumin)	neg.: 85 %, pos.: 99 %
ERY/Hb	Intact erythrocytes: 3-7 Ery/µL Hemoglobin: 5-12 Ery/µL	neg.: 99 %, pos.: 94 %

c) The values for neg. and pos. indicate the proportion of concordant negative or positive results.

Precision

Repeatability was checked for 3 test strip lots in 3 separate runs with 21 measurements per run and lot.
For details see table below.

Precision						
Repeatability with human spiked samples						
Parameter	GLU	LEU	NIT	PRO	ERY intact	ERY lysed
Target concentration[d)]	Exact agreement					
0	100.00 %	100.00 %	100.00 %	100.00 %	100.00 %	100.00 %
+	80.95 %	100.00 %	100.00 %	100.00 %	100.00 %	100.00 %
++	100.00 %	100.00 %	100.00 %	100.00 %	100.00 %	100.00 %

d) Sample concentration corresponds to the colour block
A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.
Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.
The Summary of Safety & Performance Report can be found here:
https://ec.europa.eu/tools/eudamed

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1 GP16-A3 (Urinalysis; Approved Guideline - Third Edition).
2 ECLM, European Urinalysis Guidelines. Scand J Clin Lab Invest, 2000. 60: p. 1-96.
3 Rangaiahagari A, Nyirabanzi J, Uwizeyimana JP. Comparison of urinary culture and urine dipstick nitrite test in urinary tract infection, Rwanda Medical Journal (2015);72:5-7.
4 Ronald A. The Etiology of Urinary Tract Infection: Traditional and Emerging Pathogens, Dis Mon (2003);49:71-82.
5 Susan King-Strasinger, M.S., Urinalysis and Body Fluids, 5th Edition. ISBN 978-0-8036-1697-4 (alk. paper), 2008.
6 Simerville JA., Maxted WC and Pahira J.J., Urinalysis: a comprehensive review. Am Fam Physician, 2005. 71(6): p. 1153-62.
7 Rao PV, Jones JS. How to evaluate 'dipstick hematuria': What to do before you refer. Cleveland Clinic Journal of Medicine (2008);75(3):227-233.
8 Cowart SL, Stachura ME. Glucosuria, Clinical Methods: The History, Physical, and Laboratory Examinations. 3rd edition., Boston 1990, Chapter 139.
9 Wu, A,H,B., Tietz Clinical Guide to Laboratory Tests; 4th Edition. 2006.
10 Pannala AS., et al., The effect of dietary nitrate on salivary, plasma, and urinary nitrate metabolism in humans. Free Radic Biol Med, 2003. 34(5): p. 576-84.
11 Johnson DW. Global Proteinuria Guidelines: Are We Nearly There Yet?; Clin Biochem Rev. (2011);Vol.32.
12 McPherson RA, M.R.P., HENRY'S Clinical Diagnosis and Management by Laboratory Methods, 23rd edition. ISBN 9780323295680, 2017.
13 Wollin T, Laroche B and Psooy K, Canadian guidelines for the management of asymptomatic microscopic hematuria in adults. Can Urol Assoc J, 2009. 3(1): p. 77-80.
14 Nielsen M, Qaseem A; High Value Care Task Force of the American College of Physicians. Hematuria as a marker of occult urinary tract cancer: advice for high-value care from the American College of Physicians. Ann Intern Med. 2016;164(7):488-497.
15 't Hoen LA, et al. Update of the EAU/ESPU guidelines on urinary tract infections in children. J Pediatr Urol 2021;17(2):200-207.

Symbols / Symboles / Símbolos / Simboli
Roche Diagnostics uses the following symbols and signs for this product. / Roche Diagnostics utilise les signes et les symboles suivants pour ce produit. / Para este producto, Roche Diagnostics emplea los siguientes símbolos y signos. / Per questo prodotto, Roche Diagnostics impiega i seguenti simboli:

Device for self-testing / Dispositif conçu pour l'autocontrôle / Dispositivo para autodiagnóstico / Dispositivo per test autodiagnostico

SYSTEM — Analyzers/Instruments on which reagents can be used / Analyseurs/appareils compatibles avec les réactifs / Analizadores/instrumentos adecuados para los reactivos / Analizzatori/strumenti su cui i reagenti possono essere usati

REF — Reference number / Numéro de référence / Número de referencia / Numero di catalogo

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COBAS, COMBUR-TEST and NAVIFY are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

© 2023, Roche Diagnostics

CE 0123

Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com
+800 5505 6606

05871638001
2025-01 **V 3.0**