
Cirrus™2 nebuliser Intersurgical EcoLite™ mask kits

REF	1453015, 1454015, 1453001
en	Cirrus2 nebuliser Intersurgical EcoLite mask kits
fr	Kits nébuliseur Cirrus2 avec masque Intersurgical EcoLite
de	Set mit Cirrus2 Vernebler und Intersurgical EcoLite Maske
es	Kits de nebulizador Cirrus2 con mascarilla Intersurgical EcoLite
pt	Kits de máscara Intersurgical EcoLite para nebulizador Cirrus2
it	Nebulizzatore Cirrus2, kit con maschera Intersurgical EcoLite
nl	Cirrus2 vernevelaar Intersurgical EcoLite mask sets
no	Cirrus2 nebulisator Intersurgical EcoLite maskesett
fi	Cirrus2 nebulisaattori Intersurgical EcoLite maskisetit
sv	Cirrus2 nebulisator Intersurgical EcoLite-mask set
da	Cirrus2 Intersurgical EcoLite forstøversæt med maske
el	Cirrus2 σετ νεφελοποίησης με μάσκα Intersurgical EcoLite
lt	Cirrus2 purkštuvo su Intersurgical EcoLite kauke rinkinys
pl	Zestawy nebulizatora Cirrus2 maski Intersurgical EcoLite
ru	Комплект для взрослых небулайзер Cirrus2 с маской Intersurgical EcoLite
cs	Sada Nebulizátru Cirrus2 s maskou Intersurgical EcoLite
hu	Cirrus2 porlasztó Intersurgical EcoLite maszkkészletek
sl	Cirrus2 nebulizator Intersurgical EcoLite kit z masko
lv	Cirrus2 nebulaizera Intersurgical EcoLite maskas komplekts
et	Cirrus2 nebulisaator Intersurgical EcoLite maski komplekt
bg	Cirrus2 пулверизатор Intersurgical EcoLite комплект с маска
sr	Cirrus2 nebulizator i Intersurgical EcoLite maska - set
ro	Kit nebulizare Cirrus2 cu mască EcoLite
sk	Cirrus2 Nebulizér Intersurgical EcoLite súpravy masiek
hr	Komplet raspršivača Cirrus2 i maski Intersurgical EcoLite
tr	Cirrus2 nebuliser Intersurgical EcoLite maske kitleri

1
2
3
4
Air/O2
5
6
1453015
1453001
1454015
7
1453015
1453001
1454015

en The instruction for use is not supplied individually. Only 1 copy of the instructions for use are provided per box and should therefore be kept in an accessible location for all users. More copies are available upon request.

NOTE Distribute instructions to all product locations and users. These instructions contain important information for the safe use of the product. Read these instructions for use in their entirety, including warnings and cautions before using this product. Failure to properly follow warnings, cautions and instructions could result in serious injury or death to the patient. Serious incidents to be reported as per the local regulation.

WARNING A WARNING statement provides important information about a potentially hazardous situation, which, if not avoided could result in death or serious injury.

CAUTION A CAUTION statement provides important information about a potentially hazardous situation, which if not avoided could result in minor or moderate injury.

INTENDED USE To convert a liquid drug into a mist to enable a patient to inhale and deposit the drug in the lungs.

INDICATIONS The products are designed for the therapeutic administration of nebulised medication intended for tracheobronchial deposition. The medication should be indicated for use via facemask nebulisation with a small volume jet nebuliser.

CONTRAINDICATIONS The user must follow the contraindications of the required medication. Do not use with medication that is not recommended for nebulisation via a face mask.

INTENDED USER

1453015 and 1453001 are for adult patients in prehospital, ambulance and hospital settings by appropriately qualified clinical personnel and may be used by adult homecare patients who have received instruction from appropriately qualified clinical personnel.

1454015 is for paediatric patients in prehospital, ambulance and hospital settings by appropriately qualified clinical personnel. Homecare use by paediatric patients (under 18 years) must be under the supervision of an adult who has received instruction from appropriately qualified clinical personnel.

PATIENT CATEGORY 1453015 and 1453001 are for adult patients 1454015 is for paediatric patients.

SIZE SELECTION Appropriately qualified clinical personnel should ensure that the mask fits the patient to ensure there are no leaks around the mask seal.

PRE USE CHECK Ensure all components are undamaged and are free of obstructions and foreign bodies. Do not use if packaging is damaged.

INSTRUCTIONS FOR USE

1. Wash hands and put on fresh gloves before adding medication to the nebuliser cup and administering the treatment. Use aseptic technique to avoid contamination of the medication and nebuliser cup.
2. Unscrew the nebuliser cup from the nebuliser (Fig. 1).
3. Add the medication to the nebuliser cup. The minimum fill volume when used in the vertical orientation is 2ml, the minimum fill volume when used in the horizontal is 4ml. The maximum fill volume is 10ml (Fig. 2).
4. Screw the nebuliser cup to the nebuliser, a positive click will indicate the cup is closed (Fig. 3).
5. Connect one end of the tube to the spigot connector at the bottom of the nebuliser, and connect the other end to the driving gas source -air or oxygen, as clinically appropriate. This may be a hospital pipeline system, cylinder or electrically driven compressor. When used with an air or oxygen flowmeter, the flowmeter must be pressure compensated and conform with ISO 15002. 1453001 is supplied with a 6mm flow meter nipple adaptor which can be screwed onto a 9/16 UNF Diameter Index Safety System (D.I.S.S.) flow meter to allow the oxygen tube to be attached to the D.I.S.S. flow meter. Tubing must be attached using a push/twist action (Fig. 4).
6. Attach the nebuliser to the face mask port using a push twist action (Fig. 5).
7. The adult mask elastic strap in 1453015 and 1453001 may be fitted below or above the ears. Using the below the ears position will help prevent the occurrence of pressure sores around the patient's ears. The paediatric mask elastic strap in 1454015 should be fitted above the ears only (Fig. 6).
8. Pull the elastic strap through the mask on each side and adjust for a secure, comfortable fit. Check to ensure there are no leaks around the mask, and adjust the elastic strap accordingly (Fig. 7).
9. Turn on the gas flow. The recommended gas flow required to drive the nebuliser is 6-8 litres per minute. 8 litres per minute, with the nebuliser in a vertical position, will provide optimum performance.
10. Remove the nebuliser mask kit from the patient as soon as nebulisation is completed, turn off the driving gas and remove the oxygen tubing from the driving gas source.

IN USE CHECKS

Ensure there are no leaks around the mask – adjust elastic straps if necessary.

Ensure that the tube is securely connected to the nebuliser and driving gas source.

Ensure the tubing does not cause a trip hazard.

The residual volume is the amount of medication remaining at the end of the nebulisation procedure, this may be minimised by tapping the cup when the nebuliser begins to splutter.

Check to ensure mist produced by the nebuliser is visibly entering the mask.

Monitor the patient throughout the therapy.

Ensure the mask vent slots are not occluded.

WARNINGS

1. To avoid nebulised medication entering the patient's eyes ensure there are no leaks around the mask seal.
2. Consult and follow the manufacturer's instructions for use supplied with the medication to be nebulised and observe warnings and cautions.
3. Any use of the medical device requires full understanding and observation of the instructions for use.
4. The medical device may only be used for the purpose specified as the intended use. The clinician should be aware of and consider any known pre-existing medical condition that would preclude or could affect the safe use of the product.
5. This is a single use product. Do not reuse. Reuse poses a threat to user safety in terms of cross infection, malfunction and breakage. Cleaning can lead to malfunction and breakage.
6. Reuse on a different patient poses a threat to patient safety in terms of cross infection.
7. There is a risk of fire associated with oxygen equipment. Do not use near sparks or open flames.
8. Smoking including e cigarettes/vaping during oxygen driven nebuliser therapy is dangerous and is likely to result in facial burns or death. Do not allow smoking/vaping within the same room where the oxygen source or any oxygen carrying accessories are located.
9. The recommended gas flow required to drive the nebuliser is 6-8 litres per minute.
10. Keep the tubing away from the neck to avoid airway obstruction.
11. Do not lubricate fittings, connections, tubing or other accessories to avoid the risk of fire and burns.
12. Do not use product after the expiry date.
13. The nebuliser may be used in an orientation from vertical through to horizontal only.
14. Users should follow their local and national guidelines on the use of PPE when using this product.
15. Use of this product on patients with beards, stubble or facial hair may prevent a good seal between the face and mask.
16. Cirrus 2 nebuliser technical data was obtained by nebulising 0.9% NaCl (saline) with the nebuliser in a vertical orientation. Performance characteristics can alter when nebulising different medication and by positioning the nebuliser in an alternative orientation.
17. Use only water-based lotions or salves that are oxygen-compatible prior to and during oxygen driven nebulisation therapy. Never use petroleum or oil-based lotions or salves to avoid the risk of fire and burns.

18. Oxygen makes it easier for a fire to start and spread. Do not leave the device on bed coverings or chair cushions, if the oxygen is turned on, but not in use; the oxygen will make the materials flammable. Turn the oxygen off when not in use to prevent oxygen enrichment.
19. Ensure the mask vent slots are not occluded.
20. When using the adult mask, to prevent pressure damage on the top of the ears, secure the elastic under the ears.
21. The tubing must only be used for connection between a flowmeter and an oxygen therapy or aerosol therapy device. Do not connect this tubing to a luer connector.
22. When used with an air or oxygen flowmeter, the flowmeter must be pressure compensated and conform with ISO 15002.
23. Only to be used with medication designed to be used in a small volume jet nebuliser.

CAUTION

1. For use by appropriately qualified personnel only.
2. Ensure that appropriately qualified personnel are familiar with the contents of this instruction.
3. Always perform pre-use checks.
4. FOR USA: Rx ONLY
5. Position tubing safely to avoid creating a trip hazard.
6. Cease use if the product causes any allergic reaction to the patient.
7. The user should follow local infection control guidelines and policies relating to nebuliser products.
8. Use of a driving gas flow outside the recommended minimum and maximum values will affect nebuliser performance.
9. The correct placement and positioning of the patient interface is critical to the effectiveness of the therapy.

STORAGE Recommended storage at room temperature over the indicated product shelf-life.

DISPOSAL Following use, the product must be disposed of in line with local healthcare hygiene and waste disposal protocols and regulations.

APPROPRIATE STANDARDS This product complies with relevant requirements of the following standards:

• BS EN ISO 27427 Anaesthetic and respiratory equipment. Nebulizing systems and components
• BS EN 13544-2 Respiratory therapy equipment. Tubing and connectors
• BS EN ISO 5356-1 Anaesthetic and respiratory equipment. Conical connectors. Cones and sockets
• BS EN ISO 15001 Anaesthetic and respiratory equipment. Compatibility with oxygen

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Alternative breathing patterns, including paediatric and neonatal breathing patterns, will affect nebuliser performance values. Other factors affecting nebuliser performance are the patient interface, gas flow rate, medication, nebuliser orientation and fill volume used.

Weight Intersurgical EcoLite adult mask (kg)	0.017
Weight Intersurgical EcoLite paediatric mask (kg)	0.011
Weight Oxygen tube 2.1m (kg)	0.027
Weight Oxygen Flowmeter nipple D.I.S.S. (kg)	0.002
Weight Cirrus 2 nebuliser (kg)	0.025

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Mass Median Diameter (MMD) using 8 L/min driving gas flow (µm) [1]	% particles <5 microns using 8 L/min driving gas flow [2]	Residual volume (ml) using 8 L/min driving gas flow [3]	Maximum Noise Emission (dB) [4]
3.46	66.4	1.18	<68


1. Mass Median Diameter (MMD) using 8L/min driving gas flow (µm) • Diamètre médian en masse (MMD) avec un débit de gaz moteur de 8 l/min (µm) • Mittlerer Massendurchmesser (MMD) bei einem Treibgasfluss von 8 l/min (µm) • Diámetro medio de masa (MMD) utilizando un flujo de gas de 8 l/min (µm) • Diâmetro Mediano de Massa (MMD) usando fluxo de gás de condução de 8L/min (µm) • Diametro mediano di massa (MMD) utilizzando un flusso di gas di guida di 8 l/min (µm) • Massa Mediane Diameter (MMD) bij 8 l/min gasflow (µm) • Massemediandiameter (MMD) ved bruk av 8L/min drivende gassstrøm (µm) • Massamediaanihalkaisija (MMD) käyttämällä 8 l/min ajokaasuvirtaa (µm) • Massmediandiameter (MMD) med 8L/min i flödeshastighet (µm) • Mass Median Diameter (MMD) ved hjælp af 8L/min drivende gasflow (µm) • Μέση διάμετρος μάζας (MMD) χρησιμοποιώντας ρυθμό ροής αερίου 8 L/min (µm) • Masės vidutinis skersmuo (MMD), naudojant 8 l/min. varomųjų dujų srautą (µm) • Mediana Masowa Średnicy (MMD) przy przepływie gazu napędowego 8 l/min (µm) • Массовый медианный диаметр (MMD) при расходе управляющего газа 8 л/мин (мкм) • Hmotnostní střední průměr (MMD) s použitím průtoku hnacího plynu 8 l/min (µm) • Mass Median Diameter (MMD) 8 l/perc hajtógázárammal (µm) • Srednji masni premer (MMD) z uporabo pretoka pogonskega plina 8 L/min (µm) • Masas vidējais diametrs (MMD), izmantojot 8L/min braukšanas gāzes plūsmu (µm) • Massi keskmine läbimõõt (MMD), kasutades gaasivoolu kiirust 8 l/min (µm) • Масов среден диаметър (MMD) при използване на 8L/min задвижващ газов поток (µm) • Maseni srednji prečnik (MMD) korišć enjem protoka pogonskog gasa od 8L/min (µm) • Diametrul median al masei (MMD) folosind un debit de gaz de antrenare de 8 l/min (µm) • Hmotnostný stredný priemer (MMD) s použitím prietoku hnacieho plynu 8 l/min (µm) • Srednji promjer mase (MMD) korištenjem protoka pogonskog plina od 8 L/min (µm) • 8L/dak tahrik gazı akışını (µm) kullanan Kütle Medyan Çapı (MMD)

2. % particles <5 microns using 8L/min driving gas flow • % de particules <5 microns en utilisant un débit de gaz moteur de 8 l/min • % Partikel <5 Mikrometer bei 8 l/min Treibgasfluss • % de partículas <5 micras utilizando un flujo de gas de 8 l/min • % de partículas <5 mícrons usando fluxo de gás de condução de 8L/min • % di particelle <5 micron utilizzando un flusso di gas veicolante di 8 l/min • % deeltjes <5 micron bij

gebruik van 8 l/min gasflow • % partikler <5 mikron ved brug af 8L/min drivende gasstrøm • % hiukkasia <5 mikronia 8 l/min ajokaasuvirtauksella • % partiklar <5 mikron med 8L/min i flödeshastighet • % partikler <5 mikron ved brug af 8L/min drivende gasstrøm • % σωματιδίων <5 microns χρησιμοποιώντας ρυθμό ροής αερίου 8 L/min • Dalelių procentas < 5 mikronai, naudojant 8 l /min. varomųjų dujų srautą • % cząstek <5 mikronów przy przepływie gazu napędowego 8 l/min • % частиц <5 микрон при расходе управляющего газа 8 л/мин • % částic <5 mikronů při průtoku hnacího plynu 8 l/min • 5 mikron alatti részecskék százaléka 8 l/perc hajtógázáramlással • % delcev <5 mikronov z uporabo pretoka pogonskega plina 8 L/min • Daļiņu procentuālā daļa <5 mikroni, izmantojot 8L/min braukšanas gāzes plūsmu • Osakeste % <5 mikronit, kasutades 8 l/min gaasivoolu • % частици <5 микрона при използване на 8L/min задвижващ газов поток • % čestica <5 mikrona koristeć i 8L/min protok pogonskog gasa • % particule <5 microni folosind un flux de gaz de antrenare de 8 l/min • % častíc < 5 mikrónov pri použití 8 l/min prietoku hnacieho plynu • % čestica <5 mikrona korištenjem 8 L/min pokretačkog protoka plina • % parçacıklar <5 mikron 8L/dk tahrik gazı akışı kullanılarak

3. Residual volume (ml) using 8L/min driving gas flow • Volume résiduel (ml) avec un débit de gaz moteur de 8 l/min • Restvolumen (ml) bei 8 l/min Treibgasfluss • Volumen residual (ml) con un flujo de gas de 8 l/min • Volume residual (ml) usando fluxo de gás de condução de 8L/min • Volume residuo (ml) utilizzando un flusso di gas veicolante di 8 l/min • Residuaal volume (ml) bij gebruik van 8 l/min gasflow • Restvolum (ml) ved bruk av 8L/min drivende gassstrøm • Jäännöstilavuus (ml) käyttämällä 8 l/min käyttökaasuvirtaa • Restvolym (ml) med 8L/min i flödeshastighet • Restvolumen (ml) ved brug af 8L/min drivende gasstrøm • Υπολειπόμενος όγκος (ml) χρησιμοποιώντας ρυθμό ροής αερίου 8L/min • Likutinis tūris (ml) naudojant 8 l/min. varomųjų dujų srautą • Objętość resztkowa (ml) przy przepływie gazu napędowego 8 l/min • Остаточный объем (мл) при расходе транспортировочного газа 8 л/мин • Zbytkový objem (ml) při průtoku hnacího plynu 8 l/min • Maradék térfogat (ml) 8 l/perc hajtógázárammal • Preostala prostornina (ml) z uporabo pretoka pogonskega plina 8 L/min • Atlikušais tilpums (ml), izmantojot 8L/min braukšanas gāzes plūsmu • Jääkmaht (ml), kasutades gaasivoolu 8 l/min • Остатъчен обем (ml) при използване на 8L/min задвижващ газов поток • Preostala zapremina (ml) korišć enjem protoka pogonskog gasa od 8L/min • Volumul rezidual (ml) folosind un debit de gaz de antrenare de 8 l/min • Zvyškový objem (ml) pri použití prietoku plynu 8 l/min • Preostali volumen (ml) korištenjem pokretačkog protoka plina od 8 L/min • 8L/dak tahrik gazı akışı kullanılarak kalan hacim (ml)

4. Maximum Noise Emission (dB) • Émission de bruit maximale (dB) • Maximale Geräuschemission (dB) • Emisión máxima de ruido (dB) • Emissão Máxima de Ruído (dB) • Emissione massima di rumore (dB) • Maximale geluidsemissie (dB) • Maksimal støyutslipp (dB) • Suurin melupäästö (dB) • Maximalt ljudemission (dB) • Maksimal støjemission (dB) • Μέγιστη εκπομπή θορύβου (dB) • Didžiausia triukšmo emisija (dB) • Maksymalna emisja dźwięku (dB) • Максимальный уровень шума (дБ) • Maximální emise hluku (dB) • Maximális zajkibocsátás (dB) • Največja emisija hrupa (dB) • Maksimālā trokšņa emisija (dB) • Maksimaalne müraemissioon (dB) • Максимална шумова емисия (dB) • Maksimalna emisija buke (dB) • Emisia maximă de zgomot (dB) • Maximálna hlučnosť (dB) • Maksimalna emisija buke (dB) • Maksimum Gürültü Emisyonu (dB)

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
INTERSURGICAL® COMPLETE RESPIRATORY SYSTEMS

Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
T:+44 (0)118 9656 300 info@intersurgical.com www.intersurgical.com

Ireland
T: +353 (0)1 697 8540
info@intersurgical.ie

España
T: +34 91 665 73 15
info@intersurgical-es.com

Nederland
T: +31 (0)413 243 860
info@intersurgical.nl

Sverige
T: +46 08 514 30 600
info@intersurgical.se

Россия
T: +7 495 7716809
info@intersurgical.ru

South Africa
T: +27 (0)11 444 7968
info@intersurgical.co.za

台湾
T: +886 4 2380 5430
info@intersurgical.com.tw

Canada
T: +1 905-319-6500
info@intersurgical.ca

France
T: +33 (0)1 48 76 72 30
info@intersurgical.fr

Portugal
T: +351 219 108 550
info@intersurgical.pt

België/Belgique
T: +32 5 2456991
info@intersurgical.be

Danmark
T: +45 28 69 80 00
info@intersurgical.dk

Česká Republika
T: +420 272 911 814
info@intersurgical.cz

中国
T: +86 519 69817000
info@intersurgical-cn.com

Pilipinas
T: +63 2820 4124
info@intersurgical.ph

Colombia
T: +57 1 742 5161
info@intersurgical.com

Deutschland
T: +49 (0)2241 25690
info@intersurgical.de

Italia
T: +39 0535 20836
info@intersurgical.it

Luxembourg
T: +32 5 2456991
info@intersurgical.be

Lietuva
T: +370 387 66611
info@intersurgical.lt

Türkiye
T: +90 216 468 88 28
info@intersurgical.com.tr

日本
T: 81 (0)3 6863 4359
info@intersurgical.co.jp

USA
T: +1 800 828 9633
support@intersurgicalinc.com

Australia
T: +61 (0)2 8048 3300
info@intersurgical.com.au

MADE BY INTERSURGICAL IN LITHUANIA

Distributed in the USA by Intersurgical Incorporated, 6757 Kinne Street, East Syracuse, NY 13057
T: (800) 828-9633 F: (315) 451-3696 support@intersurgicalinc.com www.intersurgical.com

2979 Issue 2 JS 05.23