
Vasofix® Braunüle®
Vasofix® Certo

B BRAUN

LLDorder 7150 - G 230706

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
5249_Vasofix_Certo_Braunüle_IFU_1123.indd 7 31.10.23 10

en Instructions for use

The Vasofix® IV Catheter with injection port is an indwelling peripheral catheter to obtain intravenous access for short-term use. It is designed to:

- Allow quick needle-free bolus injections via injection port
- Improve catheter stability and minimize catheter movement with fixation wings.

The Vasofix® IV Catheter can be used with Luer access devices which are in compliance with ISO 80369-7.

This device is intended for single use only. It is provided sterile and non-pyrogenic.

Materials Used
Vasofix® Braunüle®:
FEP, PP, PE, Silicone rubber, ABS, Chrome-Nickel steel.
Vasofix® Certo:
PUR, PP, PE, Silicone rubber, ABS, Chrome-Nickel steel.

Components are not made with natural rubber latex, PVC, or DEHP.

Indication
The Vasofix® IV Catheter is inserted into a patient´s peripheral venous system for short-term use to

- sample venous blood, or
- administer infusion solutions, blood products or pharmaceuticals (e.g. in or without carrier solutions) according to the SPC of the drug/ solution.

Patient population & Intended User
The catheter may be used for any patient population with consideration given to adequacy of vascular anatomy and appropriateness for the procedure being performed, the solution being infused and duration of therapy.

The catheter is intended to be used by qualified Healthcare Professionals who are trained on the use of peripheral IV catheters according to national regulations.

Contraindications

- The Vasofix® IV Catheter should not be used in patients with known hypersensitivity to any of the materials employed.
- The Vasofix® IV Catheter shall not be used for arterial access, subcutaneous access and high pressure power injection applications.

Precautions & Warnings

- Observe standard precautions on all patients. Aseptic technique, proper skin preparation, and continued protection of the insertion site are essential.
- When placing or maintaining any IV Catheter, standard precautions must be adhered in accordance with Centers for Disease Control and Prevention/

5249_Vasofix_Certo_Braunüle_IFU_1123.indd 8 31.10.23 10

Occupational Safety and Health Administration (CDC/ OSHA) standards for blood-borne pathogens to avoid the risk of exposure to contaminated blood.
- Use only if packaging is intact. This device is sterile unless packaging has been opened or damaged.
- Do not re-use. Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/ or limited functionality of the device may lead to injury, illness or death of the patient.
- In the case of an unsuccessful IV Catheter placement attempt, remove the needle first, then remove catheter from patient and discard both.
- Never reinsert the needle inside the catheter once the needle has been partially or completely withdrawn as it may pierce and/or severe the catheter.
- Report needlestick injuries immediately and follow established institutional protocols.
- Use of only ISO Luer slip and Luer lock connectors ensures proper functioning.
- Luer Slip connections should not be left unattended due to potential of disconnection.
- Extreme care should be taken not to damage, pierce, cut or sever the catheter. Therefore, do not bend the catheter and/or needle during insertion, advancement, or removal of the needle.
- Do not use scissors or sharp instruments at or near the insertion site.
- If blood is present, flush device according to institutional protocol.
- Monitor tightness of all connections during the whole infusion procedure.
- The catheters may be used with access pressure of up to 3 bar.
- The catheters may be used to facilitate the placement of vascular access devices such as guidewires, indwelling central venous catheters (CVC), peripherally inserted central catheters (PICC) and midline catheters into the vascular system.

Residual risks/Side effects
- General risks related to IV Catheters: infiltration, extravasation, leakage, air embolism, phlebitis, thrombophlebitis, thrombosis, catheter-related bloodstream infection, local infection, inflammation, catheter occlusion, catheter fracture, catheter embolism or clot formation.

Application
Use aseptic technique. Do not rotate catheter hub prior to insertion.
1. Gather all necessary items for insertion and stabilization procedure.
2. Select and disinfect per institutional protocol.

5249_Vasofix_Certo_Braunüle_IFU_1123.indd 9 31.10.23 10

3. Apply tourniquet.
4. Remove protective guard in a straight outward motion.
5. Anchor the vessel with gentle skin traction and access vessel. Observe blood return in flashback chamber to confirm successful vessel entry (see figure A).
6. Lower and advance the entire catheter and needle unit slightly to ensure catheter tip is in vessel (see figure B).
7. Using the injection port, advance catheter off the needle (approximately 1/8" or 3mm). You can observe blood return between the needle and catheter to confirm catheter is in the vessel (see figure C). After confirmation, continue advancing catheter off the needle and into vessel.
8. Release tourniquet.
9. Before removing the needle, occlude vein distal to catheter tip to prevent blood spillage. At the same time, stabilize the catheter hub to prevent catheter dislodgement during needle removal. Withdraw the needle straight out with a controlled and continuous motion (minimize rotation or bending of the needle) (see figure D).
10. Dispose of the needle immediately into a sharps-proof container.
11. Immediately connect the infusion line or accessory device and cover the puncture site with a sterile and transparent dressing (see figure E) per institutional protocol.
12. For bolus injections open the injection port and attach syringe to injection port with a slight twisting movement (see figure F). The port opens and closes automatically. After injection immediately close protective cap of injection port.
13. After removal, dispose catheter according to local guidelines and/ or institutional protocols.

Duration of use Duration of Use

- Duration of use is dependent on the applied therapy according to the SPC of the drug/solutions and should follow the national guidelines and/or hospital protocols.
- The insertion site should be checked frequently and at regular intervals. The catheter should be removed in event of local or systemic signs of infection.

Disposal Disposal

Disposal according to local guidelines and/ or institutional protocols.

If, during the use of this product or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

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B. Braun Melsungen AG