Product Data Sheet

1.0 Name of the product	**Leukoplast® Hospital**
2.0 Product description	
2.1 Description	Leukoplast® hospital is a medical adhesive tape that consists of a white viscose fabric coated directly with a rubber-based caoutchouc adhesive mass containing Dry Natural Rubber (DNR) and perforated by air. The device is available in different sizes
2.2 Intended purpose	Leukoplast® hospital is intended for short time use on intact skin for fixation of dressings of all types and sizes, for traction bandages or other devices. The device is non-sterile and for single application use. Intended users are health care professionals and patients. Use of the device is not restricted to specific populations.
2.3 Instructions for use	No IFU is required for this kind of product due to the low risk to user. All necessary information for safe use is available on labelling
2.4 Warnings and precautions for use	Tape contains natural rubber latex and colophony.
2.5 Contraindications	Known intolerance to latex or colophony.
2.6 Transport Precautions	not applicable
2.7 Duration of application	One day to several days depending on the skin condition.

2.9 Composition



Designation (incl. no. in drawing)	Description	% of final product
Backing (1)	*Viscose Fabric*	51
Adhesive (2)	*Natural rubber based*	49
Spool with or without snapring (3)	*Polystyrene and additive*	N/A

2.10 Product range

Size/ Color	Spools (= UoM) per Folding Box	Folding Boxes per Shipper Case	Spools per Shipper Case	Product code
9,2 m x 1,25 cm	24	10	240	01756-00
9,2 m x 2,50 cm	12	10	120	01757-00
9,2 m x 5,00 cm	6	10	60	01758-00
9,2 m x 1,25 cm	24	10	240	76820-00
9,2 m x 2,50 cm	12	10	120	76821-00
9,2 m x 5,00 cm	6	10	60	76822-00

2.12 Storage conditions

Dry, clean, and without exposure to direct sunlight

2.13 Shelf life/ Storage time

5 years

2.14 Sterilization

Not applicable

2.15 Reimbursement Information

Product is reimbursable in most countries. For more information, please reach out to your Essity contact or contact the market access department directly.

BNR	Size	PZN
01756-00	9,2 m x 1,25 cm	04593534
01757-00	9,2 m x 2,50 cm	04593540
01758-00	9,2 m x 5,00 cm	04593557
76820-00	9,2 m x 1,25 cm	08527333
76821-00	9,2 m x 2,50 cm	08527356
76822-00	9,2 m x 5,00 cm	08527362

2.16 Key quality parameters

Characteristic	Method	Target
Adhesion to steel (final product)	Essity internal	2.4 – 4.0 N/cm
Air permeability	Essity internal	25 – 250 $cm^3/cm^2/s$

3.0 Chemical substances of special concern to Essity

Essity has defined chemical substances that are of special concern and are subject to specific restrictions. A reference to the list of substances can be found in Annex A2 via the following link www.essity.com/gss

Deviations are covered in 3.1 Components.

3.1 Components

Raw materials used in product formulation

Substance	Included in formulation
natural rubber Latex	Yes
Lanolin and its derivatives	Yes
Colophony	Yes
Colophony derivatives	Yes
Bisphenol A (BPA)	No
Polyvinylchloride (PVC)	No
Methylmethacrylate	No
Butylacrylate	No
Microplastics (<5 mm)	No
Nanoparticles	No
Fragrance/parfum	No
Antibiotics	No
Triclosan	No
Chlorhexidine	No
Polyhexanide	No

Raw materials used in the packaging

Substance	Included in formulation
Polyvinylchloride (PVC)	No
Natural rubber latex	No
Recycled material	Partly in shipper case

Raw materials used in the product formulation and in the packaging may contain small amounts of the substances listed above and not marked as included in formulation as well as amounts of other substances not listed above. As trace level analysis is not performed for each batch BSN medical GmbH cannot provide proof of absence.

History

Version / Date	Page / Item	Description of Change
01 / 20.10.2023	All	Set up new document