# Discofix®

15328312
0523

en **Multidirectional stopcocks with tubing for infusion therapy**
de **Mehrwegehähne mit Schläuchen für Infusionstherapie**
bg **Многопосочни спирателни кранчета с тръбичка за инфузионна терапия**
cs **Vícesměrné uzavírací kohouty s hadičkou k infuzní léčbě**
da **Flervejsstophaner med slange til infusionsbehandling**
el **Στρόφιγγες πολλαπλών κατευθύνσεων με σωλήνα για θεραπεία έγχυσης**
es **Llaves multidireccionales con tubos para terapia de infusión**
et **Voolikuga mitmesuunalised korkkraanid infusioonraviks**
fi **Monisuuntaiset letkulliset hanat infuusiohoitoa varten**
fr **Robinets multidirectionnels avec tubulure pour perfusions**
hu **Csővezetékkel ellátott többirányú csapok infúziós terápiához**
it **Rubinetti multidirezionali con tubi per terapia infusionale**
lt **Daugiakrypčiai čiaupai su vamzdeliais infuzijų terapijai**
lv **Daudzvirzienu noslēdzošie krāni ar caurulītēm, lai veiktu infūzijas terapiju**
nl **Multidirectionele afsluitkraan met slang voor infuustherapie**
no **Flerveis-stoppekraner med rør for infusjonsterapi**
pl **Wielodrożny zawór odcinający z drenem przeznaczony do terapii infuzyjnej**
pt **Torneiras multidireccionais com tubos para terapia de infusão**
ro **Robinet de închidere multidirecțional cu tuburi pentru tratament prin perfuzie**
sk **Viaccestné ventily s hadičkami na infúznu terapiu**
sl **Večsmerni nepovratni ventili s cevkami za zdravljenje z infuzijo**
sv **Flervägsventiler med slangar för infusionsbehandling**
tr **İnfüzyon tedavisi için hatlı çok yollu musluklar**

## en Warning

- Do not use if packaging or product is damaged or contaminated or if protective caps are loose or missing. Use of damaged or contaminated device might lead to injury, illness or death of the patient.
- Do not re-use: The products are single-use, sterile and non-pyrogenic medical devices. Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and / or impairment of functional capability. Contamination and / or limited functionality of the device may lead to injury, illness or death of the patient
- Multidirectional stopcocks may tend to fissure formation when influenced as follows:
  - o getting into contact with specific infusion solutions and drugs, in particular with high PH-value or fat emulsions
  - o alcoholic disinfectants
  - o sudden changes in temperature necessitated by therapy
  - o residues of blood or thrombocytes as well as precipitates of substances physico-chemically incompatible at the junction points.

  These fissures may lead to fluid re-flux or air entry. Examine the entire infusion line regarding these hazards and observe absolutely the SPC of the respective drug manufacturer.
- Monitor all connections for tightness during the whole application.
- Do not re-attach the protective cap of the patient connector.

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B. BRAUN Vietnam Co., Ltd.
B.Braun Melsungen AG