WESTMED

Diagnostikleuchte
Diagnostic light

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HONSUN (NANTONG) Co., Ltd No. 8 Tongxing Road
Economic & Technological, Development Area 226009
Nantong City, Jiangsu Province, P. R. C.

XXXXXX

EN **Please read the instructions thoroughly and carefully before using this product.**

Application/safety: The diagnostic light may only be used by authorized specialist personnel. The instrument must be kept out of the reach of children.
Purpose: The light is used as an aid for diagnostic illumination of the area to be examined, such as the oral cavity or the skin surface.
Operation: Press the "Switch" button to switch the device on and off.
Inserting the batteries: Unscrew the lid of the diagnostic light. Now insert two AAA batteries with the positive pole towards the bulb. You can then screw the cover back in.
Changing the bulb: Unscrew the tip and remove the contact spring. Now the defective bulb can be replaced and the diagnostic light can be reassembled with the contact spring.
Warnings: Do not place the diagnostic light in liquids. Avoid direct contact with the tongue and mucous membrane during use. When examining the oral cavity, special care must be taken to ensure that the light does not fall into the oral cavity under any circumstances. There is a risk of suffocation! The trader accepts no liability for translation faults.

EC REP Shanghai International Holding Corp. GmbH (Europe)
Eiffestraße 80, 20537 Hamburg, Germany

CH REP MedNet SWISS GmbH
D4 Platz 4, 6039 Root D4 Switzerland

CE MDR (EU) 2017/745
ISO 13485:2016

REF HS-401F2 | 5810035

MD

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WESTMED® distributed by Fleischhacker GmbH & Co. KG