
i-gel® O_2

INTERSURGICAL®
COMPLETE RESPIRATORY SYSTEMS

i-gel® O_2 Resus Pack

Instructions for use

CE 1639


STERILE EO


Rx ONLY


en	i-gel® O_2 Resus Pack Instructions For Use	6 - 9
fr	i-gel® O_2 Resus Pack Mode d'emploi	10 - 13
de	i-gel® O_2 Resus-Paket Gebrauchsanweisung	14 - 17
es	Equipo de reanimación i-gel® O_2 Instrucciones de uso	18 - 21
pt	Embalagem de reanimação i-gel® O_2 Instruções de utilização	22 - 25
it	Kit di rianimazione i-gel® O_2 Istruzioni per l'uso	26 - 29
nl	i-gel® O_2 resuspakket Gebruiksaanwijzing	30 - 33
no	i-gel® O_2 resus-pakke Bruksanvisning	34 - 37
fi	i-gel® O_2 -elvytyspakkaus Käyttöohjeet	38 - 41
sv	i-gel® O_2 återupplivningskit Bruksanvisning	42 - 45
da	i-gel® O_2 Resus Pack Brugsanvisning	46 - 49
el	Σετ ανάνηψης i-gel® O_2 Οδηγίες χρήσης	50 - 53
lt	„i-gel® O_2" gaivinimo paketas Naudojimo instrukcijos	54 - 57
pl	i-gel® O_2 – zestaw do resuscytacji Instrukcja uzytkowania	58 - 61
ru	Реанимационный комплект i-gel® O_2 Инструкция по применению	62 - 65


cs	i-gel® O_2 resuscitační balíček Návod k použití	66 - 69
hu	i-gel® O_2 lélegezteto´ csomag Használati utasítás	70 - 73
sl	Pripomoček i-gel® O_2 Resus Navodila za uporabo	74 - 77
lv	Ierīces i-gel® O_2 atdzīvināšanas komplekts Lietošanas instrukcija	78 - 81
et	i-gel® O_2 elustamiskomplekt Kasutusjuhised	82 - 85
bg	Пакет i-gel® O_2 Resus Инструкции за употреба	86 - 89
sr	i-gel® O_2 pakovanje za reanimaciju Uputstvo za upotrebu	90 - 93
ro	Set de resuscitare i-gel® O_2 Instruct,iuni de utilizare	94 - 97
sk	i-gel® O_2 Resuscitačný balík Návod na použitie	98 - 101
hr	Paket za oživljavanje i-gel® O_2 Upute za uporabu	102 - 105
tr	i-gel® O_2 Resus Pack Kullanım Talimatları	106 - 109
ja	i-gel® O_2 蘇生パック 使用説明書	110 - 113
zh-s	i-gel® O_2 急救包 使用说明	114 - 117
zh-t	i-gel® O_2 復甦包 使用說明	118 - 121
ar	حزمة الإنعاش i-gel® O_2 تعليمات الاستخدام	122 - 125

Figure 1. / Figure 1. / Figure 1 / Abbildung 1. / Figura 1. / Figura 1. / Figura 1. / Figuur 1. / Figur 1. / Kuva 1. / Bild 1. / Figur 1. / Εικόνα 1. / 1 pav. / Ilustracja 1. / Рисунок 1. / Obrázek 1. / 1. ábra / Diagram 1 / 1. attēls. / Joonis 1. / Фигура 1. / Slika 1. / Figura 1. / Nákres 1. / Slika 1. / Şekil 1. / 図1 / 图 1. / 圖 1. / الشكل 1.


10

11
11c
11b
11a


12

13

14


15

16

17

18

19

20

NOTE: Distribute instructions to all product locations and users.
These instructions contain important information for the safe use of the product. Read these instructions for use in their entirety, including warnings and cautions before using this product.
WARNING: A WARNING statement provides important information about a potentially hazardous situation which, if not avoided could result in death or serious injury.
CAUTION: A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided could result in minor or moderate injury.

Product Name & Description
i-gel® O_2 Resus Pack. Product codes: 8703000, 8704000 and 8705000.

Key Components (Figure 2)
The i-gel® O_2 Resus Pack consists of the following, supplied in a sterile pack.
1. Airway support strap
2. Gastric tube
3. i-gel® O_2
4. Sachet of lubricant

Figure 2




Figure 3

Key Components
Key components (Figure 3)
1.1 Distal end of gastric channel
1.2 Epiglottic rest
1.3 Proximal end of gastric channel
1.4 Soft non-inflatable cuff
1.5 Buccal cavity stabiliser
1.6 Integral bite block
1.7 Hook ring for airway support strap
1.8 Cap for supplementary oxygen port
1.9 15mm connector
1.10 Supplementary oxygen port

Product positioning (Figure 4)
Key:
2.1 Oesophagus
2.2 Trachea
2.3 Vocal cords/folds
2.4 Thyroid cartilage
2.5 Laryngeal Inlet
2.6 Sealing mechanism. Non-inflatable cuff
2.7 Ventilatory opening
2.8 Hyoid bone
2.9 Epiglottis
2.10 Tongue
2.11 Nominal depth of insertion mark
2.12 Incisor
2.13 Nominal length of the internal pathway (refer to technical data chart below)

Figure 4

Intended Use
Securing and maintaining a patent airway. To provide a conduit for the delivery of respiratory gases to the patient.

Indications
The i-gel® O_2 is indicated for use in securing and maintaining a patent airway in routine and emergency anaesthetics of fasted patients, during spontaneous or intermittent positive pressure ventilation, during resuscitation of the unconscious patient, and as a conduit for intubation under fibreoptic guidance in a known difficult or unexpectedly difficult intubation, by personnel who are suitably trained and experienced in the use of airway management techniques and devices. The i-gel® O_2 can also be used for providing supplementary oxygen to a spontaneously breathing patient during postoperative care or patient transfer once the patient has been disconnected from the breathing circuit or passive oxygenation during resuscitation of the unconscious patient.

Contraindications
- Patients with any condition which may increase the risk of a full stomach e.g. hiatus hernia, sepsis, morbid obesity, pregnancy or a history of upper gastro-intestinal surgery etc.
- Non-fasted patients for routine and emergency anaesthetic procedures.
- Trismus, limited mouth opening, pharyngo-perilaryngeal abscess, trauma or mass.
- Do not allow peak airway pressure of ventilation to exceed 40cm H_2O.
- Do not use excessive force to insert the device or gastric tube.
- Inadequate levels of anaesthesia which may lead to coughing, bucking, excessive salivation, retching, laryngospasm or breath holding thus complicating the anaesthetic outcome.
- Do not leave the device in situ for more than 4 hours.
- Do not reuse or attempt to reprocess the i-gel® O_2.

Patient category
The product range is suitable for use with all patient populations. The product range is suitable for use with adult patients with a body weight as indicated on the individual product pack and the connector of the device. Please refer to the technical data table for size selection. Whilst size selection on a weight basis should be applicable to the majority of patients, individual anatomical variations mean the weight guidance provided should always be considered in conjunction with a clinical assessment of the patient's anatomy.

Warnings
- i-gel® O_2 must be lubricated according to the instructions for use.
- The patient should always be in the 'sniffing-the-morning air' position prior to insertion with the assistant helping to open the patient's mouth, unless head/neck movements are considered inadvisable or are contraindicated.
- The i-gel® O_2 is supplied in a protective cradle to ensure the device is retained in the correct flexion prior to use. The protective cradle also acts as a base for lubrication.
- Optimum depth of anaesthetic must be achieved prior to attempting insertion (i.e. absence of eyelash reflex, easy up and down movement of the lower jaw, no reaction to pressure applied to both angles of the mandible).
- The leading edge of the i-gel® O_2 tip must follow the curvature of the patient's hard palate upon insertion.
- If there is a failure to achieve complete insertion after utilising the standard insertion technique and a jaw thrust, deep rotation or triple maneuvre has also failed, then the device should be inserted under direct vision by laryngoscopy or one size smaller device should be used.
- After insertion, i-gel® O_2 should be taped from maxilla-to-maxilla or secured with the airway support strap provided. The i-gel® O_2 should be tied as well as taped during surgical procedures that require the patient to be in a lateral or prone position. The i-gel® may also be tied in place in circumstances where tape is considered inadequate or ineffective.
- Excessive air leak during manual ventilation is primarily due to either sub-optimal depth of anaesthesia or sub-optimal depth of i-gel® O_2 insertion.
- Particular care should be taken with patients who have an ASA or Mallampati score of III and above or who have fragile and vulnerable dental work, in accordance with recognised airway management practices and techniques.
- As with all supraglottic airways, it is important to ensure the correct size of device is used, lubrication is optimal, the device is inserted and positioned correctly and regularly checked intraoperatively in order to reduce the potential for nerve damage, tongue numbness, cyanosis and other potential complications.
- When the supplementary oxygen port of the i-gel® O_2 is being used, the patient's respiratory parameters must be continuously monitored.
- When the supplementary oxygen port of the i-gel® O_2 is being used, the proximal end of the main airway channel must be left open to air and not connected to a ventilator or Bag-Valve-Mask (BVM) resuscitation system.
- As with any airway device, when using the i-gel® O_2, it is essential to ensure there is always an adequate expiratory pathway for gas to avoid the potential for barotrauma.
- In order to reduce the risk of fire, it is essential the user fully complies with the guidance provided by the European Resuscitation Council (ERC), the American Heart Association (AHA) or other appropriate local resuscitation council regarding the safe use of oxygen during defibrillation.
- Remove dentures or removable plates from the mouth before attempting insertion of the device.
- Do not apply excessive force during insertion.
- It is not necessary to insert fingers or thumbs into the patient's mouth during the process of inserting the device.
- The i-gel® O_2 must always be separated from the protective cradle prior to insertion. The protective cradle is not an introducer and must never be inserted into the patient's mouth.
- The airway support strap should only be used on patients in a supine position. It is intended as a supplementary or short-term method of securing the i-gel® O_2 until additional securing of the device can be achieved with tape or other recognised method. The airway support strap is only suitable for use with the i-gel® O_2 and should not be used with any other device.
- The flow of oxygen through the supplementary oxygen port must not exceed 8 l/ min in order to avoid the application of peep (positive end expiratory pressure).
- No attempt should be made to use the i-gel® O_2 as a conduit for intubation without fibreoptic guidance.
- For patients undergoing Electroconvulsive Therapy (ECT) treatment, a separate bite block specifically designed to be used during ECT must always be used in conjunction with the i-gel® O_2.

Cautions
- Always perform the pre-use checks included in these instructions for use.
- The i-gel® O_2 resus pack is a single use product. Do not reuse. Reuse poses a threat to patient safety in terms of cross infection, device malfunction and breakage. Cleaning can lead to device malfunction and breakage.

Pre-use checks
- Inspect the packaging and ensure it is not damaged prior to opening.
- Inspect the device carefully, check the airway is patent and confirm there are no foreign bodies or a BOLUS of lubricant obstructing the distal opening of the airway or gastric channel.
- Carefully inspect inside the bowl of the device, ensuring surfaces are smooth and intact and also that the gastric channel is patent.
- Discard the device if the airway tube or the body of the device looks abnormal or deformed.
- Check that the 15mm connector fits the patient connection.
- Ensure the supplementary oxygen port is firmly closed with the integral cap until it is required for use.

Pre-insertion preparation
1. Always wear gloves.
2. Open the i-gel® O_2 package and take out the protective cradle containing the device **(Figure 1, picture 1)**. Remove the accessory pack containing the sachet of lubricant, gastric tube and airway support strap from the protective cradle and place to one side **(Figure 1, picture 2)**.
3. In the final minute of pre-oxygenation, remove the i-gel® O_2 from the protective cradle **(Figure 1, picture 3**) and transfer it to the palm of the same hand that is holding the protective cradle, supporting the device between the thumb and index finger **(Figure 1, picture 3)**.
4. Open the sachet of supplied lubricant and place a small bolus onto the middle of the smooth surface of the protective cradle in preparation for lubrication **(Figure 1, picture 4)**. Do not use silicone based lubricants.
5. Grasp the i-gel® O_2 along the integral bite block and lubricate the back, sides and front of the cuff with a thin layer of lubricant. This process may be repeated if lubrication is not adequate, but after lubrication has been completed, check that no BOLUS of lubricant remains in the bowl of the cuff or elsewhere on the device. Avoid touching the cuff of the device with your hands **(Figure 1, pictures 5, 6, 7 and 8)**.
6. Place the i-gel® O_2 back in to the protective cradle in preparation for insertion **(Figure 1, picture 9)**.

Recommended insertion technique
1. Remove the lubricated i-gel® O_2 from the protective cradle and grasp the device firmly along the integral bite block. Position the device so that the i-gel® O_2 cuff outlet is facing towards the chin of the patient **(Figure 1, picture 10)**.
2. The patient should be in the 'sniffing the morning air' position **(Figure 1, picture 10)** with head extended and neck flexed. The chin should be gently pressed down by an assistant before proceeding to insert the i-gel® O_2.
3. Introduce the leading soft tip into the mouth of the patient in a direction towards the hard palate.
4. Glide the device downwards and backwards along the hard palate with a continuous but gentle push until a definitive resistance is felt.
5. At this point the tip of the airway should be located into the upper oesphageal opening **(Figure 1, picture 11a)** and the cuff should be located against the laryngeal framework **(Figure 1, picture 11b)**. The incisors should be resting on the integral bite-block **(Figure 1, picture 11c)**.
6. i-gel® O_2 should be taped down from 'maxilla to maxilla' or secured with the airway support strap provided **(Figure 1, pictures 12,13,14 and 15)**.

Use of the airway support strap
An airway support strap is provided as part of the i-gel® O_2 Resus Pack. This can be used in addition to taping of the device (Recommended insertion technique – point 6), to provide additional security or as a short-term alternative to taping where the use of tape is designated difficult or unsatisfactory, such as on patients with facial hair. It can also be used as a preliminary, short-term measure to secure the device until definitive fixation can be achieved.

Recommended technique for use of the airway support strap
1. The strap should be slid under the patient's neck until the wide central band of the strap is located directly under the neck of the patient **(Figure 1, picture 13)**.
2. One end of the strap should then be lifted over the patient's face and secured to the i-gel® O_2 by placing an appropriate hole on the strap over the lug of the hook ring located at the top of the integral bite block **(Figure 1, picture 14)**.
3. The other end of the strap should then be lifted over the other side of the patient's face and secured in the same manner, ensuring there is sufficient tension to hold the i-gel® O_2 securely in place, but not an excessive tension that may cause trauma to the patient's neck or face or that may cause unwanted downward pressure of the i-gel® O_2 **(Figure 1, picture 15)**. Some adjustment of the strap may be needed to ensure optimal positioning.

Use of the gastric tube
A size 12 French gauge (FG) gastric tube is provided as part of the i-gel® O_2 Resus Pack. If required, the gastric tube can be passed down the gastric channel of the i-gel® O_2 to empty the stomach contents.

Recommended technique for inserting the gastric tube
Prior to inserting the gastric tube, ensure both the gastric tube and the gastric channel of the i-gel® O_2 have been adequately lubricated. The ideal way to achieve this is to place a small bolus of lubricant over the proximal end of the gastric channel. The gastric tube should then be inserted a short way down the gastric channel and moved up and down to prime the channel before completing insertion of the gastric tube.

Use of the supplementary oxygen port
The i-gel® O_2 consists of a standard i-gel®, with an integral supplementary oxygen port. The port is located centrally at the proximal end of the i-gel®, just above the flat end of the connector. It is supplied with an integral cap, which should remain closed when the port is not in use. The oxygen port is designed to provide supplementary oxygen to a spontaneously breathing patient during postoperative care or patient transfer once the patient has been disconnected from the breathing circuit or passive oxygenation during resuscitation of the unconscious patient. If supplementary oxygen is not required, the i-gel® O_2 can be used as a standard supraglottic airway by ensuring the supplementary oxygen port remains closed with the integral port cap.

Recommended technique for use of the supplementary oxygen port
1. With the patient disconnected from the breathing system, remove the cap from the oxygen port **(Figure 1, picture 16)**.
2. Connect one end of a standard oxygen tube to the supplementary oxygen port of the i-gel® O_2 and the other end to a suitable method of oxygen delivery **(Figure 1, pictures 17 and 18)**.
3. Set the desired flow rate on the oxygen delivery system **(Figure 1, picture 19)**.
4. When treatment has been completed, turn off the oxygen supply, remove the oxygen tubing from the i-gel® O_2 and replace the integral port cap firmly onto the spigot of the i-gel® O_2 **(Figure 1, picture 20)**.

Oxygen Flow L/Min	Approximate FiO_2 %*
4	35
6	45
8	50

*Inhaled oxygen percentage is affected by the patient's breathing pattern.

Use of the i-gel® O_2 as a conduit for intubation
In the known difficult or unanticipated difficult intubation, the i-gel® O_2 can be used for intubating the patient, by passing an endotracheal tube (ETT) through the airway channel of the device under fibreoptic guidance.

The i-gel® O_2 should always be used in accordance with recognised airway management practice for supraglottic airway devices.

In Use Checks
- The position of the device should be regularly checked intraoperatively in order to reduce the potential for nerve damage, tongue numbness, cyanosis and other potential complications.
- The gastric channel should be regularly checked for any regurgitant.
- The patency of the airway channel should be checked after any movement of the patient's head or neck.

Technical Data

Product code	Size	Patient weight (kg)	Internal volume (ml)	Nominal length of internal pathway (mm)	Maximum size of ET tube that can be inserted down the airway channel (mm)	Maximum size of gastric tube that can be passed down the gastric channel (FG)	Minimum interdental gap for insertion (mm)	Pressure drop
8703000	3	30-60	22.3	180.8	6.0	12	27.7	0.2cmH_2O @ 60L/min
8704000	4	50-90	26.5	179.2	7.0	12	30.7	0.2cmH_2O @ 60L/min
8705000	5	90+	30.4	187.2	8.0	14	32.9	0.1cmH_2O@ 60L/min

Storage Conditions
Store in a dry place out of direct sunlight.
Recommended storage at room temperature over the indicated product shelf-life.

Disposal
Following use, the product must be disposed of in accordance with local hospital infection control and waste disposal regulations.

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* Explanation of packing symbols can be found on the Intersurgical website http://www.intersurgical.com/support or in standard ISO 15223-1:2016.

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Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
T:+44 (0)118 9656 300 info@intersurgical.com www.intersurgical.com

Ireland
T: +353 (0)1 697 8540
info@intersurgical.ie

España
T: +34 91 665 73 15
info@intersurgical-es.com

Nederland
T: +31 (0)413 243 860
info@intersurgical.nl

Sverige
T: +46 08 514 30 600
info@intersurgical.se

Россия
T: +7 495 7716809
info@intersurgical.ru

South Africa
T: +27 (0)11 444 7968
info@intersurgical.co.za

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Canada
T: +1 905-319-6500
info@intersurgical.ca

France
T: +33 (0)1 48 76 72 30
info@intersurgical.fr

Portugal
T: +351 219 108 550
info@intersurgical.pt

België/Belgique
T: +32 5 2456991
info@intersurgical.be

Danmark
T: +45 28 69 80 00
info@intersurgical.dk

Česká Republika
T: +420 272 911 814
info@intersurgical.cz

中国
T: +86 519 69817000
info@intersurgical-cn.com

Pilipinas
T: +63 2820 4124
info@intersurgical.ph

Colombia
T: +57 1 742 5161
info@intersurgical.com

Deutschland
T: +49 (0)2241 25690
info@intersurgical.de

Italia
T: +39 0535 20836
info@intersurgical.it

Luxembourg
T: +32 5 2456991
info@intersurgical.be

Lietuva
T: +370 387 66611
info@intersurgical.lt

Türkiye
T: +90 216 468 88 28
info@intersurgical.com.tr

日本
T: 81 (0)3 6863 4359
info@intersurgical.co.jp

USA
T: +1 800 828 9633
support@intersurgicalinc.com

Australia
T: +61 (0)2 8048 3300
info@intersurgical.com.au

patent: http://us.intersurgical.com/support

Distributed in the USA by Intersurgical Incorporated, 6757 Kinne St., East Syracuse, NY 13057
T: (800) 828-9633 F: (315) 451-3696 support@intersurgicalinc.com www.intersurgical.com

2090 Issue 7 RO 09.21