cobas®

Roche

Urisys® 1100

Short Instructions

Publication version 1.0

Table of contents

English **3**

Español **20**

Português **39**

Urisys® 1100 Short Instructions - English

Publication information

This publication is the Short Instructions for the **Urisys® 1100**. For more information, refer to the **Urisys® 1100** Operator's Manual and the method sheets.

Electronic user documentation can be downloaded using the eLabDoc e-service on navify® Portal:
navifyportal.roche.com
Alternatively, scan the QR code.
For more information, contact your local affiliate or Roche representative.

Revision history

Publication version: 1.0
Revision date: 07-2024

Change description: For more information, refer to the **Urisys® 1100** Operator's Manual.

Incident reporting

Inform your Roche representative and your local competent authority about any serious incidents which may occur when using this product.

If you have any questions or problems, contact your Roche representative.
A list of all Roche affiliates can be found at:
www.roche.com/about/business/roche_worldwide.htm

Contact addresses

Roche Diagnostics GmbH
Sandhofer Strasse 116
68305 Mannheim
Germany
Made in Hungary

Copyright

© 2024, Roche Diagnostics GmbH.
All rights reserved.

Trademarks

COBAS, URISYS, URISYS 1100, COMBUR-TEST, and CHEMSTRIP are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Intended use

The **Urisys® 1100** is a reflectance photometer designed to read and evaluate the urine test strips Combur[10] Test® UX from Roche Diagnostics. It reads the strips under standardized conditions, saves the results to memory and outputs them via its own inbuilt printer and/or interface.
The **Urisys® 1100** is designed for In Vitro Diagnostic (IVD) use by qualified physicians and laboratory staff.
Not for self-testing.
The instrument is intended for near-patient testing.

IVD — For *in vitro* diagnostic use.

Device for near-patient testing.

CE — Complies with the provisions of the applicable EU regulations.

Device not for self-testing.

c UL us LISTED Laboratory Equipment 8C79 — Issued by Underwriters Laboratories, Inc. (UL) for Canada and the US.

Product names and descriptors

Urisys® 1100 = device
Combur[10] Test® UX = urine test strips
Chemstrip® 10 A = urine test strips
Control-Test M = calibration strips
cobas® infinity edge = service system
navify® POC Operations = data management system

Information for the Urisys® 1100

Safety instructions

⚠ General attention

To avoid serious or fatal injury, ensure that you are familiar with the system and safety information before you use the system.

- Pay particular attention to all safety precautions.
- Always follow the instructions in this publication.
- Do not use the device in a way that is not described in this publication.
- Store all publications in a safe and easily accessible place.

Safety classification

The safety precautions and important user notes are classified according to the applicable standards. Familiarize yourself with the following meanings and icons:

⚠ Safety alert

- The safety alert symbol is used to alert you to potential physical injury hazards. Obey all safety messages that follow this symbol to avoid possible damage to the system, injury, or death.

⚠ WARNING! Warning...

- ...indicates a hazardous situation that, if not avoided, could result in death or serious injury.

⚠ CAUTION! Caution...

- ...indicates a hazardous situation that, if not avoided, could result in minor or moderate injury.

NOTICE! Notice...

- ...indicates a hazardous situation which, if not avoided, may result in damage to the system.

Important information that is not safety relevant is indicated with the following icon:

Tip...

...indicates additional information on correct use or useful tips.

Safety precautions

⚠ CAUTION! Insufficient knowledge and skills

To ensure correct operation of the device:

- Only healthcare professionals may operate the device.
- Users must have received a comprehensive instruction in the operation, QC, and care of the device.

⚠ CAUTION! Non-specified consumables and accessories

Use of non-specified consumables and accessories may result in malfunction of the system and may render the warranty null and void.

- Only use consumables and accessories specified by Roche Diagnostics.

⚠ CAUTION! Personal injury and infection due to sharps, rough edges, and/or moving parts

- Good Laboratory Practice can reduce the risk of injury.
- Be aware of your laboratory environment, well-prepared, and follow the instructions for use.
- Some areas of the device may have sharps, rough edges, and/or moving parts.
 Wear personal protective equipment to minimize the risk of injury from bodily contact with such parts, especially in less accessible areas, or while cleaning the device.
- Your personal protective equipment should be appropriate to the degree and type of potential hazard, e.g. suitable lab gloves, eye protection, lab coat, and footwear.

Warning messages

⚠ WARNING! List of warning messages

Failure to observe warning messages may result in death or serious injury.

- Before using the system, read the warning messages with care.

⚠ WARNING! False test measurement

To avoid harmful outcomes resulting from false test measurements, ensure the following:

- Only use a urine sample when running a patient test.
- Regularly perform QC tests.
- Carefully follow the storage, handling, and result interpretation instructions provided by the QC material manufacturer.
- Only use consumables and accessories specified by Roche Diagnostics.
- Do not use the device outside the specified temperature range.
- Test strips and calibration strips are single-use only. Never reuse a strip.
- Do not use expired or damaged consumables.
- Carefully follow the test strip handling described in the package insert and the Operator's Manual.
- Do not use a test strip tray that is in use longer than 18 months.
- Do not touch the reference pad. Carefully handle it during cleaning procedure.
- Only use the Roche original power cable and power adapter.

⚠ WARNING! Protection against infection and pathogens

Healthcare professionals using the device to run tests must be aware that any object coming into contact with urine samples of human origin is a potential source of infection.
Users must adhere to standard precautions when handling or using the device.
All parts of this system should be considered potentially infectious and are capable of transmitting pathogens between patients, and between patients and healthcare professionals.

- Use disposable gloves.
- Wash hands thoroughly with soap and water before putting on a new pair of gloves and running the next test.
- Dispose of used test strips from patient according to your healthcare facility's infection control rules.
- Follow all local health and safety regulations.

WARNING! Avoiding electrical shock, fire, and explosions

- Only use the Roche original power cable and power adapter.
- The device must only be operated with the prescribed power adapter (Class III protection). The device is classified as Category II protection for overvoltage and Degree 2 for pollution according to IEC 664.
- Do not use loose power outlets or damaged power supplies, cables, or plugs.
- Do not short circuit the power supplies.
- Do not drop the device or power supply, and protect them against shaking and vibrations.

WARNING! Risk of injury by electronic equipment

Opening covers or removing parts of the device may expose voltage-carrying components. Connectors can be live, too.

- Do not remove any cover of the device, except where this can be achieved manually without the use of any tools.
- Never try to maintain or repair an open device which is carrying voltage.
- Do not attempt to work in any electronic compartment.

WARNING! Risk of barcode read errors

If a barcode is read incorrectly, it may lead to patient misidentification and inappropriate therapy decisions.

- When creating patient or user barcodes, always adhere to the applicable international IEC/ISO standards for the respective barcode type.
- Make sure barcode size and print quality are adequate as defined in the manual of the used barcode reader.
- Carry out a plausibility check on all data scanned into and displayed by the device.

Caution messages

CAUTION! List of caution messages

- Before operating, read the caution messages with care. Failure to observe them may result in minor or moderate injury.

CAUTION! Allergy or injury caused by reagents and other working solutions

Direct contact with reagents, detergents, cleaning solutions, disinfectants, or other working solutions may cause skin irritation or inflammation.

- Always wear disposable gloves.
- Observe the cautions given in the package inserts, or healthcare facility instructions of the reagents, cleaning solutions, and disinfectants.

CAUTION! Unauthorized system access and data loss

If the system is connected to a local area network, this network must be protected against unauthorized access. In particular, it must not be connected directly to any other network or the internet.

Users are responsible for the security of their local area network, especially in protecting it against malicious software and attacks. Failure to do so may result in data loss or may make the system unusable.

Roche recommends the following precautions:

- To protect access to all external devices, ensure to use appropriate security equipment. For more information, contact your local Roche representative.
- Make sure that the connected network is free of malicious code.
- Regularly scan attached IT systems for security vulnerabilities.

Notices

⚠ WARNING! List of notices

Failure to observe the notices may result in damage to the system.

- Before using the system, read the notices with care.

NOTICE! General care

Using cleaning solutions and disinfectants not specified in this publication may result in incorrect operation and possible failure of the device.

- Only clean the device with the recommended cleaning solutions and disinfectants.
 For more information, see the *Maintenance* section of this document.
- Make sure the device is thoroughly dried after cleaning and disinfection.

Safety messages

⚠ WARNING! Disposal of the device

Infection by a potentially biohazardous device.

- The device or its components must be treated as potentially biohazardous waste. Decontamination (i.e., a combination of processes including cleaning and disinfection) is required before reuse, recycling, returning the device to Roche, or disposal.
- Dispose of the device or its components according to the appropriate local regulations or you may return it to Roche. For more information, contact your local Roche representative.

⚠ WARNING! Access control

Unauthorized access may lead to data breaches.

- Implement physical access controls to ensure that only authorized personnel operate the system at all times.
- Assign a personal, unique user ID to each user for system access.
- Assign access rights to each user only as high as required for the tasks of the user.
- Delete user IDs from the system for users who no longer work on the system.
- Be aware of password security rules.
- Do not leave a device unlocked in publicly accessible areas.

⚠ WARNING! User authentication

The security of the system and its data depends on the password-protected access. If an unauthorized person discovers your user ID and password, they could compromise this security.

- Always enter your password unobserved.
- Do not write down your password anywhere, including in a contact form, in the address book, or in a file on the computer.
- Do not disclose your password to anyone. Roche will never ask you for your password.
- If you ever disclose your password to anyone, change it immediately afterwards.
- Contact your local Roche affiliate if you think your account has been compromised.
- Delete users who are no longer needed on the device and use the standby time function for automatic logoff on inactivity.
- Ensure that users log off the device when they have finished using it.

WARNING! Network security

Malicious software and hacker attacks may impair IT security. The laboratory is responsible for the IT security of their IT infrastructure.

- To protect and separate Roche systems from other laboratory infrastructure, the Roche-provided firewall must be used.
- Secure all devices and services used in the laboratory infrastructure against malicious software and unauthorized access.
- Secure the network environment to be resilient against traffic redirection attack and eavesdropping (man-in-the-middle attack).

WARNING! Data entry and data transfer

Writing patient-sensitive information can violate protection laws for protected health information.

- Do not write any patient-sensitive information into free-text fields.
- Data transfer using any host protocol is not encrypted; data is transferred as plaintext and readable with IT tools like network sniffers.

CAUTION! Monitor the system for suspicious activity and report suspected compromise

If you find any of the typical signs of malicious software or unauthorized access to the system (unexpected warning messages, files, or log entries like multiple failed logon attempts; significantly degraded user interface performance; seemingly random crashes of the system; automated typing of text; and so on), the following recommendations are essential:

- Physically disconnect the system from the network.
- Contact the IT staff responsible for your facility to report and verify the finding.
- Mistrust results produced while the system has been compromised.
- Contact your local Roche representative for evaluation.

CAUTION! Data loss or unavailability of the system due to malicious software or unauthorized system access

Malicious software or unauthorized system access can result in data loss or system unavailability.

To avoid infection by malicious software or the unauthorized access and misuse of the system, the following recommendations are essential:

- Make sure other computers and services on the network (for example, the LIS, archiving share, backup share, or service) are properly secured and protected against malicious software and unauthorized access.
- Restrict physical access to the system and server and all attached IT infrastructure (computer, cables, network equipment, and so on).
- Store related removable media in a safe place.
- Ensure that all device components maintaining private data are physically secure (i.e. require a physical key to be accessed or removed).
- Check all external storage devices, like CDs, DVDs, or USB flash drives with an antivirus software (on another computer) to ensure that they are free of malicious software, before using them on the system.
- Implement security and authentication measures for removable media (i.e. USB flash drives, CD-ROM/DVD, tapes, etc.).
- Use state-of-the-art security mechanisms (e.g. WPA2 EAP, if supported) to protect Wi-Fi connections, if you use a Wi-Fi adapter.
- Allow only defined users to read and write encryption keys (IIS App config/Windows Server hardening).
- Allow only internet access to trusted web sites and web services.
- Enable Data Execution Prevention for all applications (Windows Server).
- Make sure to use secure channels to download software updates of the system.

CAUTION! Unauthorized access to the system configuration

Unauthorized access to the system configuration (administration interface, configuration store, query packages) can lead to system failure, corrupted data, and unauthorized disclosure of protected health information.

- Make sure that your Wi-Fi connections are properly secured.
 - If supported, use EAP security for Wi-Fi connections.

CAUTION! Unauthorized access to confidential data

Unencrypted communication can lead to unauthorized access to confidential data.

- Use encrypted communication channels from the application.
- Contact the IT staff responsible for your facility when a security warning related to the digital certificate is displayed.

CAUTION! Sensitive data - data tampering in transmission

Unprotected Wi-Fi connections between devices and the application can lead to unauthorized data access and tampering.

- Make sure that your Wi-Fi connections are properly secured, if you use a Wi-Fi adapter.
- If supported, use EAP security for Wi-Fi connections.

CAUTION! Data loss and disaster recovery

Unexpected disasters, such as hard disk failures, system software failures, or natural disasters such as fires or floods can occur.

- The IT staff responsible for your facility must develop a disaster recovery plan to protect the entire system.
- Ensure that backups are made regularly (ideally every day) and stored securely.

NOTICE! Data transmission is not secured

Unsecured communication can compromise the integrity and confidentiality of the transmitted data. It is the customer's responsibility to ensure that the network infrastructure provides appropriate protection for the transmitted data.

- Such protection can be achieved by system level tunneling protocols or other means of network encryption.

NOTICE! Data is not encrypted

- Data in the database is not generally encrypted. It is the customer's responsibility to ensure that access to the data storage and backup files is permitted only to authorized entities.

WARNING! Secure data storage

Unauthorized access to data backups and archive files can violate data protection laws.

- Any data backup or data archive that has been exported from the device must be physically stored in a secured location.
- Ensure only authorized persons may access the secure data storage. This includes the data transfer to remote storage locations and disaster recovery.
- Data backups must not be taken from the secure data storage. Do not take external storage devices outside the laboratory environment.

WARNING! Cybersecurity and privacy awareness

Insufficiently informed employees can endanger security.

- Perform regular cybersecurity and privacy awareness trainings for staff handling personal data. Instruct staff how to handle data in a compliant way and according to the privacy principles as mandated by customer regulations.
- Check your device for suspicious activity and report any suspected compromise to your local Roche representative immediately.
- Update to the latest software versions provided by Roche as soon as possible.
- Carefully use external storage devices.
 Do not connect to the system any external storage device that you use on public or home computers.
 Failure to do so may result in data loss and render the device unusable.

WARNING! Use of storage devices

Wrong handling of external storage devices may result in data loss or system malfunction.

- Only insert or remove an external storage device after the data export is finished.
- At any one time only one external storage device can be in use. Before inserting an external storage device into a USB port, check that no other external storage device is connected.

WARNING! Data backup

Data may get lost due to memory failures.

- Back up your data (measurement results and system parameters) at regular intervals.
- Do the backup periodically to store relevant data on an external storage device.
- Make a backup copy if you have changed any system parameters.

WARNING! Non-approved third-party software

Installation of any third-party software that is not approved by Roche Diagnostics may result in incorrect behavior by the system.

- Do not copy or install any software or software patches on the system unless it is part of the system software or your local Roche representative advises it.

WARNING! Read the Operator's Manual carefully before installation to ensure proper operation of the device from the start.

NOTICE! If your device has been exposed to considerable changes in temperature and/or humidity, wait for it to acclimatize sufficiently before using it.

NOTICE! Upon receipt, check the contents of the packaging for completeness. If the contents have suffered any damage in transport, contact your local Roche representative.

CAUTION! To ensure your readings are always accurate, do not set up the device near devices that create high-frequency fields, as they may interfere and produce false results. Such devices include, for instance, walkie-talkies, mobile telephones, microwave ovens, and diathermic equipment.

NOTICE! To ensure that urinalysis is carried out correctly, read the package insert included in the test strip package.

WARNING! When measuring in Fast mode, allow the strips to incubate for approximately 45 seconds before starting the measurement. If the incubation time is too short, false-low, or false-negative results may be measured for some parameters. If the incubation time outside the device is too long, false-high results may be measured for some parameters.

WARNING! Thermal printing paper is sensitive to light and subject to fade with prolonged exposure to bright light. Patient reports are, thus, best kept in a place away from direct light and high temperature. For anticipated storage beyond 5 years it is recommended that reports are photocopied for storage.

CAUTION! Potentially incorrect results

Regular calibration is necessary to ensure the quality of the results obtained. Roche Diagnostics cannot warrant the correctness of results if the device is not calibrated regularly.

- Perform calibrations regularly.

CAUTION! Incorrect results due to biased calibration results

Soiling on calibration strips may influence the calibration results and therefore impair the validity of the test results.

- Do not touch the pads of the calibration strips.
- Do not bend the calibration strip during operation.
- Do not re-use calibration strips. Always use a new calibration strip for each calibration, including repeats. Dispose of the strip after calibration.

CAUTION! Potentially incorrect results

Regular QC measurements are necessary to ensure the quality of the results obtained. Roche Diagnostics cannot warrant the correctness of results if QC measurements are not performed regularly.

- Perform QC measurements regularly.

NOTICE! To avoid false positive results, do not pipette or drop the QC materials onto the test pads.

NOTICE! Risk of unauthorized access

Resetting the device to default settings will preserve existing user and password information, but will allow an automatic logon (without authentication) with administrator rights after the reset.

- Make sure to set the **ID setup** to the desired options right after the reset, if the device should be further used.

WARNING! Potentially incorrect results

Changing the range limits to other values than the preset defaults may have an influence on the correctness of results and is not warranted by Roche Diagnostics.

- You are responsible for validating the consistency of results after making changes to the range limits.

WARNING! Protection against infection

Healthcare professionals using the device must be aware that any object coming into contact with urine samples of human origin is a potential source of infection. Users must adhere to standard precautions when handling or using the device.

- Use disposable gloves.
- Use cleaning and disinfection materials of recommended quality only.
- Dispose of used cleaning materials according to your facility's guidelines.
- Follow all health and safety regulations in force locally.

NOTICE! Damage to the device by using unsuitable cleaning materials

Using unsuitable cleaning materials may result in incorrect operation and possible failure of the system.

- Clean and disinfect the device with recommended materials only. Make sure that the device is thoroughly dried after cleaning and disinfecting.
- Do not use any abrasive cleaning substances or pointed implements to clean the device.
- Do not spray anything onto the device and do not immerse it in liquid.

NOTICE! Damage to the device by liquids reaching internal parts

Liquids entering internal parts may lead to device malfunction.

- Only wipe the housing, never use sprays.
- Make sure that no liquid enters internal components.
- Do not immerse the device or any parts of it.
- Always ensure that the device is thoroughly dried after cleaning and disinfecting.

Environmental conditions

Temperature range	Operation: +15 °C to +32 °C Storage: -5 °C to +45 °C Transport: -25 °C to +70 °C
Relative humidity	Operation: 20% to 80% Storage and transport: up to 85%
Maximum altitude	2000 m
Extraneous light	Without influence on measured values: 20 000 lux daylight, 1000 lux artificial light
Location of operation	The device is designed for indoor use only. Do not expose the device to direct sunlight or another source of direct light (for example a spot lamp).

Requirements

Power connection	• Consumption: Operating max. 18 W; Standby 1.3 W • The device must only be operated with the prescribed power adapter (Class III protection). The device is classified as Category II protection for overvoltage and Degree 2 for pollution according to IEC 664.

Operation

Roche recommends to regularly check for software updates as described in the Operator's Manual.

Installation

- Unpack the system and place it on a stable, level surface.
- Insert the test strip tray and the printer paper roll.
- Connect the power cable to the power adapter and to a readily accessible mains connector.
- Connect the power adapter to the USB-C port at the rear of the device.

Powering on the device

- At the back of the device, press the power button. The device is powered on.
- If you install the device for the first time, or reset the device settings, choose the language of the device.
- The device performs a self-check. On completion of the self-check, the test strip tray is transported to the start position and the retaining bar opens.
- A printout containing the **Self-check OK** info and the time stamp is generated. Depending on your defined user identification option, the **User logon** or the **Home** screen is displayed.
- If you install the device for the first time, or reset the device settings, perform a calibration as requested by the device. To perform calibration, you must firstly enter the calibration strip lot information on the device.

Waking the device from standby mode

Depending on your **Power management** settings, the device switches to standby mode after the defined time of inactivity.

- Tap on the screen. The device wakes from standby mode. Depending on the device's configuration, the **User logon** screen or the previously active screen is displayed.

Logging on

Depending on your defined user identification option, you are either automatically logged on or you have to enter your credentials.

For more information on how to define new users and user identification options, see the **User management** section in the Operator's Manual.

First log on
If you log on for the first time, use the credentials (user ID and initial password) you were given by your administrator. Change the initial password to your personal password.

To log on using auto logon

- Power on the device or wake it from standby mode. You are automatically logged on. The **Home** screen or your previously active screen is displayed.

To log on manually

- On the **User logon** screen, enter your user ID or choose it from the user list and tap the ✓ button. As an operator, if you enabled the **No password required** option, you are logged on. Administrators must always enter a password, except when **Auto logon (Administrator)** is enabled.
- Enter your password and tap the ✓ button.
 When entering your password, the characters are displayed as asterisks (*). By choosing the 👁 button, you can make your password visible to, e.g., check for correctness. To hide your password again, choose the 👁 button again.
 The **Home** screen or your previously active screen is displayed.
- Make sure to remember your password. If you cannot remember your password, it can be reset by a user with administrator rights.

Entering information about consumables

The device stores and uses information about consumables used for patient tests, QC, and calibration. You have to enter this information on the device before performing any measurement.


Calibration

The device is calibrated before leaving the factory. When installed, it must be recalibrated with the specified calibration strip before the first samples are read, and thereafter every 7 days.

Performing a calibration

To perform a calibration

- When the notification about a required calibration (**No valid user calibration**) is displayed, to directly start the calibration process, tap the ✓ button.
 The **User calibration** screen is displayed.
 You can postpone a calibration by enabling the **Postpone** option before tapping the ✓ button. After postponing a calibration, the previously active screen is displayed again and a calibration will not be started. Postponing a calibration is valid for the current login session.
- To start a calibration without the notification being displayed, do the following: On the **Home** screen, tap the **Calibration** button.
 The **User calibration** screen is displayed.
- To select the required calibration lot, tap the yellow-highlighted **Lot No.** entry field.
- To select a calibration strip lot, tap the calibration lot entry, or tap the ▲ and ▼ buttons (if more than three lots are available).
- To confirm the selection and return to the **User calibration** screen, tap the ✓ button.
- On the **User calibration** screen, to start the calibration process, tap the ✓ button.
 A callout instructs you to place the calibration strip.
- Take a calibration strip from the calibration strip container.
 Be careful not to touch the test pads and do not allow them to come into contact with urine.

- Place the calibration strip, with the test pads facing upwards, on the test strip tray so that its leading edge is held by the clip (**A**) at the front end of the insertion slot. The retaining bar (**B**) must be open.
 Do not bend the calibration strip.
- To confirm the information screen and to proceed with the calibration process, tap the ✓ button.
 After warming up time, the test strip tray advances slightly. The retaining bar closes, and the grey reference pad on the test strip tray and the calibration pads are read.
 The test strip tray is then transported back to the start position, and the retaining bar opens.
- When the calibration process is completed, a corresponding message is displayed and the calibration report together with the date and time is printed.
 To display the calibration report, tap the ✓ button.
 If the calibration fails and an error message is displayed, see: Error messages and troubleshooting
- Remove and dispose of the calibration strip.
 Use each calibration strip only once.


Performing a measurement

The device can be operated in 2 different modes:


- **Normal mode**
- **Fast mode**. You operate in **Fast mode** if you run a serial measurement for which test strips dipped in urine incubate outside the device. In this mode, the throughput is increased.
 For more information, refer to the Operator's Manual.

To perform a measurement


- On the **Home** screen, tap the **Patient test** button.
- On the **Patient test** screen, tap the yellow-highlighted **Pat. ID** entry field to enter a patient ID.
 If you work with automatic patient IDs, the entry field will be automatically filled, but can be edited manually.
- CAUTION! Do not write any patient-sensitive information in patient IDs.
 Enter the patient ID using the virtual keyboard or scan it using the barcode reader.
 To save and return to the **Patient test** screen, tap the ✓ button.
 The maximum length of a patient ID is 20 characters. You will not be able to enter any more than this. If you scan a patient ID longer than 20 characters, it is truncated to the first 20 characters.
- On the **Patient test** screen, tap the yellow-highlighted **Lot No.** entry field to choose the strip lot you want to use.
- To choose the strip lot you want to use, tap the ✓ button.
 If more than 3 lots are available, tap the ▲ and ▼ buttons to scroll.
- Tap the ✓ button.
 You return to the **Patient test** test screen.
- NOTICE! To avoid incorrect readings due to a discoloration of the test strip pads, the test strip container must be closed immediately after removal of a test strip, using the original desiccant-filled stopper.
 Take a test strip from the test strip container.
- Dip the test strip briefly (1 s) in the urine sample.


- CAUTION! Do not bend the test strip.
 Draw the long edge of the test strip along the rim of the sample container to remove excess urine.

- Place the test strip, with the test pads facing upwards, on the test strip tray so that its leading edge is held by the clip (**A**) at the front end of the insertion slot. The retaining bar (**B**) must be open.
 About 2 mm of the test strip must be held under the clip (**A**). Do not bend the test strip. It is important that the test strip is correctly positioned and ready to be read within 5-10 seconds of dipping the test strip.

- On the **Patient test** screen, to start the measurement process, tap the ✓ button.
 After warming up, the test strip tray advances slightly, the retaining bar (**B**) closes, and the grey reference pad on the test strip tray is read.

- CAUTION! If the test strip is not correctly located in the middle of the test strip tray, move it gently to the side until it is properly aligned. Be careful not to move the test strip tray. Ensure that the retaining bar (**B**) is locked into place and that the test strip is in the correct position.
- 55 seconds after the measurement process has been started, the first test pad is measured, followed by the others. After that, the test strip tray returns to the start position and the retaining bar opens. The result will be displayed and, depending on your printing settings, printed, and, depending on your connection settings, transferred.
- Remove and dispose of the test strip.
- Wipe off any urine residues from the test strip tray with a lint-free cloth.
- On the **Patient test** screen displaying the test results, to return to the **Home** screen, tap the ⌂ button.

Result handling

- After each measurement, the test results are displayed (and stored) together with the following information: sequence number (**Seq. No.**), patient ID (**Pat. ID**), date and time, strip type used for measurement, lot No. of used strip type, and flags (e.g., for pending calibration or positive test results).
- Test results which diverge from negative or normal values are: Flagged with an asterisk before the parameter concerned, highlighted red.
- A red asterisk displayed in the line below the heading but before the first parameter of the test results list indicates that the reading was taken with changed range limits.
- For test results that were measured after you have postponed a due calibration, the note **(Calibration pending)** is displayed.
- For test results that were measured in **Fast mode**, the note (**Fast mode**) is displayed.


Patient report

- The patient report contains the test results and is printed together with the following information: device name, serial number of the device (**Ser. No.:**), user ID (depending on your logon settings), date and time, sequence number (**Seq. No.:**), patient ID (**Pat:**).
- Test results which diverge from negative or normal values are flagged with an asterisk before the parameter concerned.
- An asterisk printed in the line below the heading but before the first parameter of the test results list indicates that the reading was taken with changed range limits.
- If the measurement was performed with a pending calibration, **Repeat calibration** is printed below the test results.

Performing a QC measurement

To perform a QC measurement

- On the **Home** screen, tap the **QC** button.
- On the **QC** screen, tap the **QC measurement** button.
- On the **QC measurement** screen, to select the QC lot and the level you want to use, tap the QC lot/level entry, or tap the ▲ and ▼ buttons (if more than 4 entries are available).
- To confirm the selection and to proceed with the test, tap the ✓ button.
- On the **QC measurement** screen tap the yellow-highlighted **Lot No.** entry field for the test strip lot, if you want to change the test strip lot. To confirm the selection and to proceed with the test, tap the ✓ button.
- NOTICE! To avoid incorrect readings due to a discoloration of the test strip pads, the test strip container must be closed immediately after removal of a test strip, using the original desiccant-filled stopper.
 Take a test strip from the test strip container.
- Dip the test strip briefly (1 s) in the QC sample.
- CAUTION! Do not bend the test strip.
 Draw the long edge of the test strip along the rim of the sample container to remove excess control material.

- Place the test strip, with the test pads facing upwards, on the test strip tray so that its leading edge is held by the clip (**A**) at the front end of the insertion slot. The retaining bar (**B**) must be open. About 2 mm of the test strip must be held under the clip (**A**). Do not bend the test strip. It is important that the test strip is correctly positioned and ready to be read within 5-10 seconds of dipping the test strip.
- To start the measurement process, tap the ✓ button.
 The measurement process is performed as described in section: Performing a measurement
- After the measurement process, the result is displayed and printed. In case of an out-of-range measurement, **QC fail** is displayed. The parameter concerned is flagged with an asterisk and highlighted red.
- Remove and dispose of the test strip.
- On the **Last QC result** screen, you have the following options: To display the QC results of the last QC measurements for the individual levels, tap the **L1**, **L2** or **L3** button. To transfer the displayed results, tap the button.
 To print the displayed results, tap the button.
 To return to the **Home** screen, tap the button.
 To return to the **QC** screen and perform another QC test (e.g., with a different level), tap the ← button.
 To scroll through the results, use the ▲ and ▼ buttons on the right side.
- Wipe off any urine residues from the test strip tray with a lint-free cloth.

Logging off

To log off

- On the **Home** screen, tap the **Log off** button. The user is logged off and the **User logon** screen is displayed.

If **Auto logoff timeout** is enabled, the logged on user is automatically logged off after the configured time of inactivity.

Powering off the device

To power off the device

- On the **Home** screen, tap the ⏻ button. The device is powered off.

Configuration

Refer to the Operator's Manual for the following topics:

- Device settings
- Measurement settings
- ID setup
- User management
- Strip lots and calibration lots
- Connecting the device to a wired or wireless network, or a USB device such as an external keyboard, or a barcode reader

Maintenance

You must regularly clean and disinfect the exterior parts and surfaces of the device as well as the test strip tray.

Only use recommended cleaning solutions and disinfectants to clean and disinfect the device. Cleaning solutions: commercial (neutral) cleaning solution, tap water. Disinfectants: Isopropanol 70%, Ethanol 70%, Bacillol AF.

For more information on maintenance, refer to the Operator's Manual.

Troubleshooting

For more information on error messages and troubleshooting, refer to the Operator's Manual.

Available accessories and consumables

	Product name	Description	Catalogue number
Package content	**Urisys® 1100**	• Device • Power adapter, USB-C • Power cable (Euro plug to IEC 60320 C7) • Test strip tray, type C • Printer paper • Short Instructions	REF 03617548XXX
Accessories	Spare reference pads (5 pieces)	Spare reference pads	REF 11907131001
	Test strip tray, type C	Test strip tray	REF 03666735001
	Datalogic QD2590	2D hand-held barcode reader with Datalogic USB cable	REF 09831053001
Consumables	Thermal paper (20 rolls)	Printer paper	REF 11907131001
	Combur[10] Test® UX (100 test strips)	Test strips	REF 11544373XXX
	Chemstrip® 10 A (100 test strips)	Test strips	REF 11379208119
	Control-Test M (50 calibration strips)	Calibration strips	REF 11379194263
Spare parts	Power cable	Power cable (Euro plug to IEC 60320 C7)	REF 04340639001
	Power supply unit USB-C/5V/3.6A	Power adapter with USB-C output without power cable	REF 09963138001

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Copyright © 2024, Roche Diagnostics GmbH. Todos os direitos reservados.
COBAS, URISYS, URISYS 1100, COMBUR-TEST and CHEMSTRIP are trademarks of Roche.