DATA SHEET WESTMED® WESTMED® Pressure infusion cuff

WESTMED



Order number
5810000

Description

Our WESTMED® pressure infusion cuff is used for pressure-supported infusion and transfusion. The cuff is suitable for 500 ml infusion liquids and is equipped with a precision manometer and pump ball. The manometer has a scale from 0 to 300 mmHg. The pressure can be regulated using the valve on the pump ball.

Technische Informationen	
Reference number	HS-20GO
Dimensions	33 x 17 cm
Measuring range	0 - 300 mmHg
Scalability	2 mmHg
Display precision	± 3 mmHg (temperature: 18 °C – 35 °C) ± 6 mmHg (temperature: 36 °C – 50 °C)
Air leak	< 4 mmHg/min
Use temperature	0 - 50 °C
Storage temperature	-20°C bis 70°C at a relative humidity of 10 - 85% (non-condensing)
Shipping conditions	Protected against light and water
Shelf Life	10,000 measurements or 5 years from production date
Packaging unit	1 piece per package 50 packages per outer carton 16 outer cartons per pallet

©Fleischhacker GmbH & Co. KG · An der Silberkuhle 18 · 58239 Schwerte · T +49 (0) 2304 931 0 · info@fleischhacker.biz
WESTMED® is a registered word/figurative mark of Fleischhacker GmbH & Co. KG

WESTMED

General information	
Safety instructions	• Do not use if the blood pressure monitor is damaged or dirty. • Never expose the cuffs to intensive solar radiation • Never touch the cuff fabric or parts with a sharp instrument, since this could cause damage • Always deflate cuff completely before storage. • Do not dismantle manometer under any circumstance.
Cleaning an desinfection	Remove the bladder first and wipe the velcro, bladder and tubes with a damp cloth. The cuff can be washed with soap and cold water like all the other cuffs, but you must rinse the cuffs with clear water afterwards and let them air dry. Wipe off the manometer and bulb with a damp cloth. Sterilization is not necessary, since the parts do not come into direct contact with the patient‘s body.

Regulatory information	
Classification	Class I medical device According to Annex VIII, Regulation (EU) 2017/745

Verpackungssymbole	
MD	Medical device
REF	Reference number/article number
EC REP	Authorized representative Europe
CH REP	Authorized representative Switzerland
	Distributor
	Importer
	Manufacturer
	Date of production
	Observe the operating instructions
	Keep dry
	Protect against heat

©Fleischhacker GmbH & Co. KG · An der Silberkuhle 18 · 58239 Schwerte · T +49 (0) 2304 931 0 · info@fleischhacker.biz
WESTMED® is a registered word/figurative mark of Fleischhacker GmbH & Co. KG

CE	Die CE-Kennzeichnung steht für „Conformité Européenne“, was „Europäische Konformität“ bedeutet. Sie symbolisiert die Konformität des Produktes mit den geltenden Anforderungen, welche die Europäische Gemeinschaft an den Hersteller stellt.

HONSUN (NANTONG) Co.,
Ltd No. 8 Tongxing Road Economic & Technological
Development Area, 226009 Nantong City, Jiangsu Province, P. R. China

EC REP Shanhai International Holding Corp. GmbH (Europe)
Eiffestraße 80, Hamburg, Germany

CH REP MedNet SWISS GmbH
D4 Platz 4, 6039 Root D4, Switzerland

©Fleischhacker GmbH & Co. KG · An der Silberkuhle 18 · 58239 Schwerte · T +49 (0) 2304 931 0 · info@fleischhacker.biz
WESTMED® is a registered word/figurative mark of Fleischhacker GmbH & Co. KG