
Product information

Reference Number	BD Vacutainer® Safety-Lok™ Blood Collection Sets	Box/Case
367281	With Luer, 21 gauge 3/4 in, 12 in Tube	50/200
367287	With Luer, 21 gauge 3/4 in, 7 in Tube	50/200
367283	With Luer, 23 gauge 3/4 in, 12 in Tube	50/200
367292	With Luer, 23 gauge 3/4 in, 7 in Tube	50/200
367285	With Luer, 25 gauge 3/4 in, 12 in Tube	50/200
367294	With Luer, 25 gauge 3/4 in, 7 in Tube	50/200
367296	Without Luer, 21 gauge 3/4 in, 12 in Tube	50/200
367297	Without Luer, 23 gauge 3/4 in, 12 in Tube	50/200
367298	Without Luer, 25 gauge 3/4 in, 12 in Tube	50/200
	BD Vacutainer® Holder	
364879	BD Vacutainer® Holder (Yellow)	25/200
368872	BD Vacutainer® Pronto™ Quick Release Holder	25/200
	BD Vacutainer® Adapters	
364902	Luer-Lok™ Access Device	50/200
367290	Multiple Sample Luer Adapter	100/1000
	BD Vacutainer® Tourniquet	
367205	Stretch Latex-free Tourniquet	25/500







Disclaimer: The information provided herein is not meant to be used, nor should it be used, to diagnose or treat any medical condition. All content, including text, graphics, images and information etc., Contained in or available through this literature is for general information purposes only. For diagnosis or treatment of any medical condition, please consult your physician/doctor. Becton Dickinson India Private Limited and or its affiliates and its employees are not liable for any damages/claims to any person in any manner whatsoever.

BD Life Sciences – Specimen Management
Becton Dickinson India Private Limited
5th & 6th Floor, Signature Tower -B,
South City - 1, Gurgaon-122001
Haryana, India
Tel: 91-124-4124300, Fax: 91-124-2383224/5/6

bd.com

©2021 BD. BD, the BD Logo and Vacutainer are trademarks of Becton, Dickinson and Company. Printed in India

BD Vacutainer®
Safety-Lok™ Blood Collection Set

● Convenience ● Compatible ● Commitment

Convenience

because practicing safety should be simple

- The BD Vacutainer® Safety-Lok™ Shield enhances safety — yet requires no change in technique.
- The slim, easy-to-handle design enhances control and flexibility.
- Large wing size permits improved positioning during use.
- Minimized coiling of tubing facilitates handling.
- Blood flashback is seen clearly through translucent shield.
- Convenient, integrated multiple sample luer adapter allows specimen collection into a "closed system" — further enhancing safety in use.

Compatible

because your focus should be on your patient

- Two tubing sizes and three needle gauges to adapt to sampling situations
- BD BACTEC™ Blood Culture Bottles are compatible with BD Vacutainer® holders allowing blood sampling from the same venipuncture


Commitment

because your product should be reliable

- The Safety-Lok™ Shield reduces needlestick risk throughout the hospital.
- After use, the needle can be easily and completely covered by the protective shield and securely locked in place.
- Improved safety from activation of the safety shield through final disposal in an approved sharps collector—may help reduce long-term costs of needlestick injuries.


Safely shielded with the

BD Vacutainer® Safety-Lok Blood Collection Set




Points to Practise

STEP 1




Apply light pressure to site using three fingers, as shown.

STEP 2




Withdraw blood collection set by grasping the translucent yellow safety shield grip area with the thumb and index finger.

STEP 3




With opposite hand, grasp tubing between thumb and index finger.

STEP 4



Push the yellow shield forward until the safety shield is locked in place.

STEP 5



Dispose of the BD Vacutainer® Safety-Lok™ Blood Collection Set and holder in the nearest sharps container.